LIQUID HANDLING TECHNOLOGY

Mixer HC

User Manual
Bedienungsanleitung
Manuel d'utilisation
Istruzioni per l'uso

EC Declaration of Conformity

USA Scientific certifies that

Mixer HC · Cat. No. 8012-0000 including Thermoblocks

has been manufactured and inspected in accordance
with the following European standards:

2014/35/EU (Low Voltage Directive)
EN 61010-1:2010 (Safety requirements for electrical equipment for control
and laboratory use)
EN 61010-2-010:2014 (Heating of materials)
EN 61010-2-051:2015 (Mixing and stirring)
2014/30/EU (EMC Directive)
EN 55011:2011 (Radio frequency disturbance)
EN 61326-1:2013 (Lab Equipm. – General EMC Requirements)
2011/65/EU (RoHS Directive)
2012/19/EU (WEEE Directive)
EC 1907/2006 (REACH Directive)
94/62/EC (Packaging and Packaging Waste Directive)
EN 50419:2006-06 (Marking of electrical and electronic equipment)
EN50581: 2013
UL 61010-1:2012, 3rd Ed. (Electrical equipment for control, and laboratory use)
CAN/CSA-C22.2 No. 61010-1-2012 (Electr. Safety Lab. Equipm.)

This product is manufactured in adherence to
ISO 9001 and ISO 13485 standards.

USA Scientific, Inc. · O Box 3565 · Ocala, FL 34478-3565
www.usascientific.com

User Manual	3
Bedienungsanleitung	31
Manuel d'utilisation	59
Istruzioni per l'uso	87
Ordering information / Bestellinformationen / References de commande / Informazioni per l'ordine	114
1 Operating instructions	5
1.1 Using this manual	5
1.2 Warning signs and hazard icons	5
1.2.1 Warning signs	5
1.3 Symbols used	5
1.4 Abbreviations used	6
1.5 Glossary	6
2 Product description	7
2.1 Main illustration	7
2.2 Delivery package	7
2.3 Features	8
2.3.1 Mixer HC	8
2.3.2 ISO Racks (optional accessory)	9
3 Safety	10
3.1 Intended use	10
3.2 User profile	10
3.3 Note on product liability	10
3.4 Warnings for intended use	10
3.5 Warning and instruction signs on the device	13
4 Installation	14
4.1 Preparing installation	14
4.2 Selecting location	14
4.3 Installing the unit	14
5 Operation	15
5.1 Overview of operating controls	15
5.2 Inserting exchangeable thermoblocks	16
5.2.1 Exchangeable thermoblocks for micro test tubes	17
5.2.2 MTP exchangeable thermoblock	17
5.2.3 Adapter plate for 0.2 ml tubes or 96 Well PCR plates	18
5.3 Inserting tubes and plates	19

5.4 Time counting mode	19
5.5 Temperature control and mixing	20
5.5.1 Basic temperature control	20
5.5.2 Operation without temperature control	20
5.5.3 Mixing	20
5.5.4 Short Mix	21
5.5.5 Interrupting operation	21
5.5.6 Interval Mix	21
5.5.7 Programmable temperature control	22
5.5.8 Save the values	23
6 Troubleshooting	23
6.1 General errors	23
6.2 Error messages	24
7 Maintenance	25
7.1 Cleaning	25
7.2 Disinfection / Decontamination	26
7.3 Decontamination before shipping	26
8 Transport, storage and disposal	27
8.1 Transport	27
8.2 Storage	27
8.3 Disposal	27
9 Technical data	28
9.1 Power supply	28
9.2 Ambient conditions	28
9.3 Weight / dimensions	28
9.4 Application parameters	28
9.5 Additional data	29
9.5.1 Interfaces	29
9.5.2 Time interval	29

1 OPERATING INSTRUCTIONS

1.1 Using this manual

- › Before using the device for the first time, please read the user manual.
- › Please view this user manual as part of the product and keep it somewhere easily accessible.
- › If the device is transferred to a third party, include this user manual.

1.2 Warning signs and hazard icons

1.2.1 Warning Signs

Depiction	Meaning
	DANGER Risk of electric shock with potential for severe injury or death as a consequence.
	DANGER Risk of explosion with potential for severe injury or death as a consequence.
	WARNING Warning of potential injury or health risk.
	CAUTION Refers to risk of damage to property.
	Refers to particularly useful information and tips.

1.3 Symbols used

Depiction	Meaning
›	You are requested to perform an action.
1. 2.	Perform these actions in the sequence described.
•	List.
Text	Terms from the device display.

1.4 Abbreviations used

DWP	Deepwell plate
MTP	Micro test plate
PCR	Polymerase Chain Reaction
Rpm	Revolutions per minute

1.5 Glossary


Deepwell plate Plate with 48, 96 or 384 wells with a larger volume than microplates. Suitable for the preparation, mixing, centrifuging, transporting and storing of solid and liquid samples.

Microplate Plates with 24, 48, 96 or 384 wells for the preparation, mixing, centrifuging, transporting and storing of solid and liquid samples.

Semi-skirted PCR plate PCR plate with surrounding half-edge.

Skirted PCR plate PCR plate with a surrounding full-edge.

Non-skirted PCR plate PCR plate with no surrounding edge.

Well Cavity. Microplate, PCR plate or Deepwell plate tube.

2 PRODUCT DESCRIPTION

2.1 Main illustration


Fig. 1 Mixer HC

1 Display	5 Mains connection socket
2 Thermosocket	6 Fuse connector
3 Operating controls	7 Mains switch
4 RS 232 interface	

2.2 Delivery package

Quantity	Description
1	Mixer HC without exchangeable thermoblock
1	Mains cable with plug
1	Hexagon socket key
1	User Manual
1	Quick-Start Guide

2.3 Features

2.3.1 Mixer HC

The robust, space-saving design of the Mixer HC combined with the chemically-resistant plastic housing and the keypad, makes this product another professional companion in your daily routine.

The Mixer HC offers many options for temperature control and mixing of liquids in the microliter and milliliter ranges:

- In exchangeable thermoblocks, liquids can be simultaneously mixed and temperature-controlled up to 99 °C
 - in 0.2, 0.5, 1.5, 2.0 ml microcentrifuge tubes
 - in 5.0 ml preparation tubes
 - in 5.0 ml test tubes
 - in conical centrifuge tubes (15 ml and 50 ml)
 - in 1.5 and 2.0 ml Cryogenic tubes,
 - as well as PCR plates, or MTP or DWP with any bottom shape style
- Cooling down the samples to 13 °C below room temperature is achieved by Peltier cooling.

Two special functions extend the area of application of the Mixer HC

- Interval mixing:
Continuous alternation between mixing phase and standstill where mixing frequency and duration of the mixing and pause functions can be freely selected.
- Program mode:
In addition to a normal mixing/temperature control run, two combined, consecutive mixing/temperature control runs can be freely programmed. When the first run is completed, the second run will switch on automatically so that two different treatments of the same liquid can be carried out.

2.3.2 ISO Racks (optional accessory)

With the ISO Racks, a total of 24 microcentrifuge tubes can be placed on an exchangeable thermoblock. This means that with the Mixer HC and the ISO Rack, microcentrifuge tubes can be handled on the lab bench at all the temperatures relevant for everyday lab use.

In addition, the tube holder serves as a working aid during the filling of microcentrifuge tubes as well as the preparation, sorting, distribution, transporting, storing and freezing of samples in microcentrifuge tubes

ISO Rack

The ISO Rack holds 24 microcentrifuge tubes (0.5 ml or 1.5/2.0 ml) for transfer of the tubes to the exchangeable thermoblocks. The lockable and stackable snap-together ISO Racks are also particularly suited to the storage of microcentrifuge tubes in refrigerators or freezers. The ISO Racks can be autoclaved. The spacing of the tubes in the ISO Rack is determined by the use of multi-channel pipettes.


Fig. 2 Removing the ISO Rack lid:

- Hold the lid on the corners and remove it in an upwards motion. If you hold the lid in the middle, it cannot be removed.

For further details please contact your local distributor.

3 SAFETY

3.1 Intended use

The Mixer HC is intended exclusively for indoor use and is for the temperature control and mixing of aqueous solutions in closed, micro tubes (0.2 to 5.0 ml) and plates.

Use only accessories recommended by the manufacturer.

For further details contact your local distributor.

3.2 User profile

This unit must only be used by skilled personnel with the appropriate training.

Before using the unit, read the user manual carefully and familiarise yourself with the unit's mode of operation.

3.3 Note on product liability

In the following cases the protection provided in the unit may be impaired. The liability for the function of the unit passes to the operator if:

- The unit is not used in accordance with the user manual.
- The unit is used outside the range of application described herein.
- The unit is used with accessories or consumables (eg, tubes and plates) which are not recommended by the manufacturer.
- The unit is maintained or repaired by persons not authorised by the manufacturer.
- Unauthorised modifications have been made to the unit.

3.4 Warnings for intended use

Please read the user manual and note the following general safety instructions before using the Mixer HC.

Risk of explosion!

- Do not operate the unit in rooms where work is being carried out with explosive substances.
 - Do not use this unit to process any explosive or highly reactive substances.
 - Do not use this unit to process any substances which could create an explosive atmosphere.
-

Danger! Electric shock from damage to unit/power cable.

- › Do not switch on the unit if either it or the power cable is damaged in anyway.
- › Do not use power cables, which were not supplied by manufacturer.
- › Only use the unit if it has been properly installed or repaired.

Danger! Lethal voltages inside the unit.

- › Ensure that the housing is always closed and undamaged so that no parts inside the unit can be accidentally touched.
- › Do not remove the unit's housing.
- › Do not allow the unit to be opened by anyone except service personnel who have been specifically authorised by your local distributor.

Danger! Electric shock as a result of penetration of liquid.

- › Switch off the unit and disconnect from the power supply before cleaning or disinfecting.
- › Do not allow any liquids to enter the inside of the housing.
- › Do not use spray disinfection.
- › Only reconnect the unit to the power supply when it is completely dry.

Danger! Damage to health from biologically or chemically hazardous substances.

- › Hazardous chemicals cause burns and other health hazards.
- › Observe the material safety data sheets for the biological and chemical substances used.
Wear personal protective equipment (PPE) at all times when working with biological or chemical substances.
- › Follow the instructions for cleaning and decontamination, and ensure that hygiene safety standards are maintained.

Risk when handling toxic or radioactively-marked liquids or pathogenic germs.

- › Follow national regulations governing the handling of these substances.
- › For complete instructions regarding the handling of germs or biological material of risk group II or higher, please refer to the „Laboratory Biosafety Manual“ (Source: World Health Organization, current edition of the Laboratory Biosafety Manual).

Warning! Damage due to incorrect power supply.

- › Only connect the unit to power sources that match the electrical specifications on the unit's nameplate.
- › Use only sockets with protective earth.

Warning! Risk of burns from hot surfaces.

The exchangeable thermoblock and the thermal base can be very hot after heating and cause burns.

- Allow heated exchangeable thermoblocks and the thermal base to cool down completely before removing the exchangeable thermoblock.

Warning! Risk of fire caused by flammable liquids.

The heating function / hot surfaces of the Mixer HC might ignite flammable liquids.

- Do not use the Mixer HC for mixing flammable liquids.

Warning! Risk of injury through chemical or mechanical damage to exchangeable thermoblocks.

- Check the exchangeable thermoblocks regularly.
- Do not use exchangeable thermoblocks that show signs of corrosion or mechanical damage.

Warning! Poor safety due to incorrect accessories and spare parts.

The use of accessories and spare parts other than those recommended by the manufacturer may impair the safety, function and precision of the unit. The manufacturer accepts no warranty or liability for damage caused by non-recommended accessories and spare parts or incorrect use.

- Use only accessories recommended by the manufacturer and original spare parts. For details, please contact your local distributor.

Warning! Risk to health from contaminated device.

- Perform decontamination before storing or dispatching the unit and/or its accessories.

Warning! Contamination due to opening tube lids.

In the following cases, the lids of tubes can spring open and the sample material can escape:

- high vapor pressure of the content,
- improperly sealed lid,
- damaged sealing lip.
- Always check that all tubes have been sealed tightly before use.

Caution! Strong vibration.

Items near the unit may be moved by vibrations on the work surface when the unit is using high speeds and may fall off the bench.

- Do not place easily movable items near the mixer or secure them adequately.
-

Caution! Damage to the display due to mechanical pressure.

- Do not exert mechanical pressure on the display.
-

Caution! Damage to electronic components from condensation.

- After moving the unit from a cooler environment (eg, cool room), wait at least an hour before connecting it to the mains power supply.
-

Caution when using aggressive chemicals.

Aggressive chemicals may damage both the unit and its accessories.

- Do not use any aggressive chemicals, such as strong or weak bases, strong acids, acetone, formaldehyde, halogenated hydrocarbons or phenol.
 - If the unit becomes contaminated with aggressive chemicals, clean it immediately with a neutral cleaning agent.
-

Caution! Lack of safety due to missing user manual.

- When passing on the unit, always enclose the user manual.
 - The current version of the operating manual and the safety instructions can also be obtained from your local distributor.
-

3.5 Warning and instruction signs on the device

Depiction	Meaning	Location
	Warning! Risk of burns from hot surfaces. The exchangeable thermoblock and the thermal base can be very hot after heating and cause burns. ➤ Allow heated exchangeable thermoblocks and the thermal base to cool down completely before removing the exchangeable thermoblock.	Top of device

4 INSTALLATION

4.1 Preparing installation

Retain all packaging for subsequent safe transport or storage.

- › Check the completeness of delivery according to 2.2 Delivery Package.
- › Check all parts for any damage during transit.

4.2 Selecting location

Select a location for the Mixer HC in accordance with the following criteria:

- Mains power connection (230 V/120 V) as per device identification plate, located on the rear of the device.
- At least 10 cm away from other equipment and walls.
- Solid bench with stable, even work surface.

4.3 Installing the unit

Warning! Damage due to incorrect power supply.

- › Only connect the unit to power sources that match the electrical specifications on the device's nameplate.
 - › Use only sockets with protective earth.
-

1. Place the Mixer HC on a suitable work surface so that the ventilation slots on the base of the unit are not blocked.
2. Connect the power cable supplied to the mains connection socket (5, Fig. 1) of the Mixer HC and the power supply.
3. Switch on the Mixer HC using the mains switch (7, Fig. 1).
4. Fit the exchangeable thermoblock to the thermosocket (2, Fig. 1).

5 OPERATION

5.1 Overview of operating controls

Before using the Mixer HC for the first time, familiarise yourself with the operating controls and the display.


Fig. 3 Operating controls

1 Setting the temperature	5 Starting interval mix
2 Setting the mixing frequency	6 Starting or stopping a run
3 Setting the run time	7 Starting Short Mix
4 Programming mixing/temperature control runs	


Fig. 4: Display

1 Nominal temperature	4 Display interval mixing time
2 Actual temperature	5 Display interval mixing function
3 Mixing frequency	6 Total time

After switching on, the display shows the values of the last run.

5.2 Inserting exchangeable thermoblocks

Warning! Risk of injury through chemical or mechanical damage to exchangeable thermoblocks.

- › Check the exchangeable thermoblocks regularly.
- › Do not use exchangeable thermoblocks that show signs of corrosion or mechanical damage.

Warning! Injury or property damage from improperly attached exchangeable thermoblocks.

If the exchangeable thermoblocks are not attached properly, they can become detached from the device or break.

- › Tighten the screws of the exchangeable thermoblocks so they are hand-tight.

Warning! Contamination due to opening tube lids.

In the following cases, the lids of tubes can spring open and the sample material can escape:

- › high vapor pressure of the content,
 - › improperly sealed lid,
 - › damaged sealing lip.
 - › Always check that all tubes have been sealed tightly before use.
-


The following exchangeable thermoblocks can be set on the Mixer HC:

Thermoblock for	24 x 0.5 ml microcentrifuge tubes
Thermoblock for	24 x 1.5 ml microcentrifuge tubes
Thermoblock for	24 x 2.0 ml microcentrifuge tubes
Thermoblock for	24 x micro tubes, diameter 12 mm
Thermoblock for	8 x 5.0 ml preparation tubes
Thermoblock for	8 x 15 ml conical centrifuge tubes
Thermoblock for	4 x 50 ml conical centrifuge tubes
Thermoblock for	24 x 1.5/2.0 ml cryogenic tubes
Thermoblock for	MTP / DWP with lid
Adapter plate	96 x 0.2 ml for PCR plates and tubes for MTP / DWP thermoblock

To insert the exchangeable thermoblocks, use the accessories supplied (eg. screw driver, screws, flat washers and lock washers).

5.2.1 Exchangeable thermoblocks for micro test tubes (microcentrifuge / cryogenic / micro tubes)

1. Place the exchangeable thermoblock (2) on the thermosocket with both flat edges and the writing facing the front.
2. Tighten the captive screws (1).


5.2.2 MTP exchangeable thermoblock

Caution! Microplates melt at extremely high temperatures.

Only heat microplates made of polystyrene to a maximum 70 °C.

1. Insert the short screws in the MTP exchangeable thermoblock with the flat washers.
2. Secure the screws from below using the lock washers.
The lock washers prevent the screws from falling out of the bores if the exchangeable thermoblock is not fixed on the unit.
3. Place the MTP thermoblock on the device in such a way that the spring plate is to the front left.
4. Tighten the screws.

5.2.3 Adapter plate for 0.2 ml tubes or 96-well PCR plates

1. Remove the lock washers (4) from the short screws in the exchangeable thermoblock.
2. Remove the short screws.
3. Connect the adapter plate (2) to the MTP thermoblock (3) using the long screws (1).
4. Secure the screws from below using the lock washers.
The lock washers prevent the screws from falling out of the bores if the exchangeable thermoblock is not fixed on the unit.
5. Place the adapter plate on the device with the MTP thermoblock.
6. Tighten the screws.

5.3 Inserting tubes and plates

Caution! Microplates melt at extremely high temperatures.

Only heat microplates made of polystyrene to a maximum 70 °C.

1. Select a suitable exchangeable thermoblock and mount it (refer to page 16).
2. Place the tubes or plate on the thermoblock.

For optimum temperature control, transfer the micro test tubes to the exchangeable thermoblocks using ISO Racks. The ISO Rack reduces the heat exchange of the thermoblock with the environment. If necessary close the ISO Rack lid.

3. If using the microplate thermoblock: put on the lid.

5.4 Time counting mode

The time counting function enables temperature control and/or mixing for a limited period of time. The unit will hold a constant temperature even after the set time has elapsed.

There are two ways to measure the time of the temperature control processes:

- **Time control:**

Time counting begins immediately when "START/STOP" is pressed. The red indicator lamp in the key lights up.

- **Temp control:**

Time counting only starts when the nominal temperature value has been reached.

While the unit heats up to reach the nominal temperature, the colon of the time display will flash. This signals that the unit is active. When the nominal temperature has been reached and time counting starts, the red indicator lamp of the "START/STOP" key lights up. To terminate the process early, press the "START/STOP" again.

1. Hold the "PROGR." key while switching on the unit. Hold the key until the display either shows **time control** or **temp control**.
2. Using the "TIME" arrow keys, set time counting.

22/24 °C	1100rpm
1h:00m	□□off

On the display, the time will appear on the bottom left. The time can be set from 1 minute to 99 hours 59 minutes. Pressing the key once will alter the setting by one minute. During the run the remaining run time will be displayed in minutes. The last minute is counted down in seconds.

5.5 Temperature control and mixing

After switching on, the values used the last time the unit was used before switching off will be shown on the display. The actual temperature flashes. The unit is now in basic temperature control mode.

The Mixer HC can cool to 13 °C below room temperature.

5.5.1 Basic temperature control

- Set the temperature with the "TEMP" arrow keys or press "START/STOP", if you do not want to change the nominal temperature.

22/24 °C	1100rpm
1h:00m	□□off

The device starts the temperature control. On the display, the nominal temperature is displayed on the left, the actual temperature on the right. The temperature between 1°C and 99°C is then adjusted immediately. When the nominal value is reached, only this value will be displayed.

The temperature can be adjusted gradually by pressing the key briefly, or adjusted continuously by pressing and holding the key.

You can also change the nominal temperature while the unit is running. The run will be continued with the re-set nominal temperature.

5.5.2 Operation without temperature control

- Press the "TEMP" arrow key until off appears on the display.
The unit remains at room temperature and does not heat. You can now use the unit as a mixer without temperature control.

5.5.3 Mixing

The mixing frequency (-- rpm) is displayed on the right of the display next to the temperature.

22/24 °C	1100rpm
1h:00m	□□off

1. Set the mixing frequency with the „MIX“ arrow keys.
You can set the mixing frequency to 0 or between 300 and 1500 rpm. Pressing the key once will adjust the frequency by 50 rpm.
Depending on the software version, the unit will automatically recognise the mounted exchangeable thermoblock. The mixing frequency specified for the exchangeable thermoblock cannot be increased in this case.

- Press the "START/STOP" key to start the mixing process.
During the run, the green indicator lamp in the "START/STOP" key lights up.
- To finish the mixing process, press "START/STOP" again.
After setting the required parameters, the unit will start mixing and/or the temperature control after pressing the "START/STOP" key. The temperature control may have also started by re-setting the nominal value for the temperature control.
When the pre-set time has elapsed or the "START/STOP" key is pressed again, mixing is completed. The temperature is, however, kept constant.

5.5.4 Short Mix

- Keep the "SHORT MIX" key pressed
The unit mixes with the set frequency until the "SHORT MIX" key is released. During the first minute, time counting starts in second intervals and then changes to minute intervals.

5.5.5 Interrupting operation

- Keep the "START/STOP" key pressed for more than two seconds.
Mixing and time counting are interrupted. You can now carry out manual operations (eg. adding reagents, replacing tubes).
- To continue the run, press the "START/STOP" key again.

5.5.6 Interval Mix

The "INTERVAL MIX" function enables a continuous alternation of mixing activity and mixing pause. The device carries out the mixing procedure, interrupted by pauses, until the set total time has elapsed. The times for mixing and the mixing pause (from 3 seconds to 99 minutes 59 seconds) can be set independently

- Press the "INTERVAL MIX" key for two seconds or more.

22/20 °C	1100rpm
1h:00m	□ 0m:5s

The indicator lamp of the key lights up in red. On the display the □ sign appears. The interval function is now activated.

- Using the "TIME" arrow keys, set the total time for the interval function.
- Briefly press the "INTERVAL MIX".
The □ sign flashes on the display.

22/20 °C	1100rpm
1h:30m	◻ 0m:10s

- Using the "TIME" arrow keys, set the mixing phase.

- Briefly press the "INTERVAL MIX" key.

The  sign flashes on the display.

- Using the "TIME" arrow keys, set the pause phase.
- Press the "START/STOP" key to start the interval function.
The unit begins to alternate between mixing and mixing pause until the total time has elapsed.
- To terminate the function before the time is up, hold the "INTERVAL MIX" key for more than 2 seconds.

The indicator lamp goes out and the display shows **Off** .

5.5.7 Programmable temperature control


The function "Programmable temperature control" enables a program run which consists of up to two steps. The second step will automatically start after the first step has been completed.

- Hold the "PROGR." key for more than two seconds.

The red indicator lamp in the key lights up. On the display the sign appears to the left of the temperature display.

- Using the arrow keys, set the parameters. The interval mixing function can be used in addition.
- Briefly press the "PROGR." key.
On the display the  sign appears to the left of the temperature display. You are now in the second stage of the program.

- Using the arrow keys, set the parameters. The interval mixing function can be used in addition.
- Press the "START/STOP" key to start the program.
The program step which is currently visible on the display is now started. If it is the first step, the second step is automatically started after the first step has been completed.

To end the program before the run is complete:

1. Press the "START/STOP" key again.
2. To return to the mode of basic temperature control, press the „PROGR.“ key.
In the „Temp control“ mode the  signs will appear for the individual programming steps.

5.5.8 Save the values

The set values will be saved even after the unit has been switched off.

6 TROUBLESHOOTING

6.1 General errors

If the proposed measures to remedy the fault repeatedly do not deliver the desired result, contact your local distributor.

Symptom / message	Cause	Remedy
Display remains dark.	No mains connection. Defective fuse.	➤ Check the mains connection and the power supply. ➤ Switch the device on. ➤ Change the fuse.
Nominal temperature is not reached.	At nominal values over 4 °C above room temperature.	➤ Set up the device in a cooler environment.
The unit does not mix or control the temperature.	Various causes are possible.	➤ Contact your local distributor.

6.2 Error messages

Symptom / message	Cause	Remedy
MOTOR ERR	Motor is not turning.	1. Switch unit off and wait for 10 seconds. 2. Switch the device on.
SPEED ERR	Motor is turning at over-speed.	1. Switch unit off and wait for 10 seconds. 2. Switch the device on.
TOO HOT!	Thermoblock hotter than 110 °C.	1. Switch unit off and wait for 10 seconds. 2. Switch the device on.
T-SENSOR ERR	Temperature sensor error.	1. Switch unit off and wait for 10 seconds. 2. Switch the device on.
TEMP ERR!	Nominal temperature is not reached.	1. Switch unit off and wait for 10 seconds. 2. Switch the device on.
EEChecksumERR RAM ERR EEPROM ERR	Memory error.	1. Switch unit off and wait for 10 seconds. 2. Switch the device on.
HEAVY LOAD	The motor does not reach nominal speed.	➤ Check if the movement of the exchangeable thermoblock is obstructed (eg. by jammed object).

7 MAINTENANCE

7.1 Cleaning

Regularly clean the housing of the Mixer HC and the exchangeable thermoblocks. For further information and detailed instructions for cleaning and disinfection contact your local distributor.

Danger! Electric shock as a result of penetration of liquid.

- › Switch off the unit and disconnect it from the power supply before starting to clean or disinfect.
 - › Do not allow any liquids to enter the inside of the housing.
 - › Do not use spray disinfection.
 - › Only reconnect the unit to the power supply once it is completely dry.
-

Caution when using aggressive chemicals.

Aggressive chemicals may damage both the unit and its accessories.

- › Do not use any aggressive chemicals on the unit and accessories such as strong or weak bases, strong acids, acetone, formaldehyde, halogenated hydrocarbons or phenol.
 - › If the unit becomes contaminated with aggressive chemicals, clean it immediately with a neutral cleaning agent.
-

Caution! Corrosion from aggressive cleaning agents and disinfectants.

- › Do not use corrosive cleaning agents, aggressive solvents or abrasive polishes.
-

Auxiliary aids

- Lint-free cloth.
- Mild, soap-based household cleaning agent.
- Distilled water.

Cleaning the Mixer HC

- Switch off the Mixer HC and isolate from the power supply.
- Clean all of the outer parts of the Mixer HC with a mild, soap-based solution, using a lint-free cloth.
- Wipe off the soap solution with distilled water.
- Dry all cleaned parts.

7.2 Disinfection / Decontamination

Danger! Electric shock as a result of penetration of liquid.

- Switch off the unit and disconnect it from the power supply before starting cleaning or disinfecting.
 - Do not allow any liquids to enter the inside of the housing.
 - Do not use spray disinfection.
 - Only reconnect the unit to the power supply once it is completely dry.
-

Auxiliary equipment

- Lint-free cloth.
- Disinfectant.

Disinfecting the Mixer HC

- Switch off the Mixer HC and isolate from the power supply.
- Allow the unit to cool down.
- Disinfect the unit.
- Select a disinfection method which complies with the legal requirements and regulations applicable to your range of application.
- Wipe the surfaces with the lint-free cloth and disinfectant.

7.3 Decontamination before shipping

If you are shipping the unit to the authorised Technical Service for repairs or to your local distributor for disposal, please note the following:

Warning! Risk to health from contaminated unit

1. Follow the instructions for the decontamination supplied by your local distributor.
 2. Decontaminate all the parts you want to dispatch.
 3. Enclose the fully-completed decontamination certificate for returned goods (incl. the serial number of the unit) with the dispatch.
-

8 TRANSPORT, STORAGE AND DISPOSAL

8.1 Transport

- Only transport the unit in the original packaging.

	Air temperature	Rel. humidity	Air pressure
General transportation	-20 to 60 °C	10 to 95 %	30 to 106 kPa
Air freight	-20 to 55 °C	10 to 95 %	30 to 106 kPa

8.2 Storage

	Air temperature	Rel. humidity	Air pressure
In transport packaging	-25 to 55 °C	10 to 95 %	70 to 106 kPa
Without transport packaging	-5 to 45 °C	10 to 95 %	70 to 106 kPa

8.3 Disposal

In the event of disposing of the product, please observe the applicable legal regulations.

Information on the disposal of electrical and electronic devices in the European Community:

The disposal of electrical devices is regulated within the European Community by national regulations based on EU Directive 2012/19/EU pertaining to waste electrical and electronic equipment (WEEE).

In accordance with this, any devices delivered after 13/08/2005 on a business-to-business basis, which includes this product, may no longer be disposed of in household waste. To document this they have been marked with the following identification:

Because disposal regulations may differ from country to country within the EU, please contact your local distributor for information.

9 TECHNICAL DATA

9.1 Power supply

Mains connection:	100 to 240 V ± 10 %, 50 to 60 Hz
Fuses:	See nameplate or label on the rear of device
Power consumption:	90 W

9.2 Ambient conditions

Environment:	For indoor use only
Ambient temperature:	4 to 35 °C
Relative humidity:	70 % max.

9.3 Weight / dimensions

Dimensions (W x D x H):	220 x 250 x 125 mm
Weight basic device:	3.2 kg

9.4 Application parameters

Temperature control range:	
Exchangeable thermoblocks for microcentrifuge tubes	13 °C below room temperature to 99 °C
Exchangeable thermoblocks for micro test plates	10 °C below room temperature to 99 °C
Temperature accuracy for exchangeable thermoblocks:	
for 1.5 ml microcentrifuge tubes:	
at nominal values between 20 °C and 45 °C	± 0.5 °C
at nominal values between <20 °C and >45 °C	± 2.0 °C
Exchangeable thermoblocks for micro test plates	
at nominal values <70 °C	± 2.0 °C
at nominal values between 70 °C and 99 °C	± 5.0 °C

Heating rate:	approx. 5 °C/min NOTE: temperature changes in filled tubes may be slower.
Cooling rate:	2 to 3 °C/min
Cooling rate from room temperature to 13 °C below room temp:	0.5 to 1.0 °C/min
Mixing frequencies for thermoblocks:	
0.5 ml microcentrifuge tubes	300 to 1500 1/min
1.5 ml & 2.0 ml microcentrifuge tubes	300 to 1400 1/min
Cryo tubes	300 to 1400 1/min
Micro test plates	300 to 1400 1/min 300 to 800 1/min when loading >200 g
15 ml & 50 ml conical centrifuge tubes	300 to 750 1/min
5.0 ml preparation tubes	300 to 750 1/min

9.5 Additional data

9.5.1 Interfaces

EDP connection:	RS-232, Sub-D9 male
-----------------	---------------------

Only connect devices to interfaces that meet the IEC 950/EN 60950 (UL 1950) standards.

9.5.2 Time interval

Programmable time interval:	1 min to 99:59 hours, infinitely adjustable
-----------------------------	---

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