Microcentrifuge 24

24-place, air-cooled Microcentrifuge

User Manual
Mode d'emploi
Manual de instrucciones

Certificate of compliance

USA Scientific certifies that
Representative samples of

Microcentrifuge 24 - Cat. No. 8030-0000 / 8030-0001 incl. rotors

Have been investigated by UL in accordance with the Standard(s) indicated
on this Certificate according to the current UL requirements.

Standard(s) for Safety:

UL 61010-1, 3rd Edition, May 11, 2012,
Revised April 29, 2016,
CAN/CSA-C22.2 No. 61010-1-12, 3rd Edition,
Revision dated April 29, 2016,
IEC 61010-1:2010 (Third Edition)

Additional Standards:

IEC 61010-2-020:2016;
IEC 61010-2-101:2015

Additional Information:

See the UL Online Certifications Directory at
www.ul.com/database for additional information.

Only those products bearing the UL Certification Mark
should be considered as being covered by
UL's Certification and Follow-Up Service.

USA Scientific, Inc.
PO Box 3565, Ocala, FL 34478-3565
www.usascientific.com

Public Health
England

Public Health England
National Infection Service
Porton Down
Salisbury
Wiltshire
SP4 0JG

Certificate of Containment Testing

Containment Testing of Rotor standard-24-at in a USA Scientific Microcentrifuge 24

Report No. 18/008 B

Report Prepared For: USA Scientific

Issue Date: 26 June 2019 re-issue

Test Summary

Rotor standard-24-at was containment tested in a USA Scientific Microcentrifuge 24, using Annex AA of IEC 61010-2-020:2016 (3rd Ed.). The sealed rotor was shown to contain a spill.

Report Written By

Name: Ms Anna Moy

Title: Biosafety Scientist

Report Authorised By

Name: Mrs Sara Speight

Title: Senior Biosafety Scientist

Please be aware that the use of the Royal Coat of Arms is highly restricted and cannot be copied. Please do not put the PHE logo on your website or use our name to endorse your products. Any reference to PHE needs to be approved by us before it can be used.

User Manual	5–33
Mode d'emploi	35–63
Manual de instrucciones	65–93

Certificates	3–4
--------------	-----

1 Operating instructions	7
1.1 Using this manual	7
1.2 Danger symbols and danger levels	7
1.2.1 Danger symbols	7
1.2.2 Danger levels	7
1.3 Symbols used	7
1.4 Abbreviations used	7
2 Safety	8
2.1 Intended use	8
2.2 User profile	8
2.3 Information on product liability	8
2.4 Application limits	8
2.4.1 Declaration concerning the ATEX directive (2014/34/EU)	8
2.5 Warnings for intended use	9
2.5.1 Personal injury or damage to device	9
2.5.2 Incorrect handling of the centrifuge	10
2.5.3 Incorrect handling of the rotors	10
2.5.4 Extreme strain on the centrifugation tubes	11
2.6 Safety instructions on device and accessories	11
3 Product description	12
3.1 Product overview	12
3.2 Delivery package	12
3.3 Features	13
3.4 Name plate	13
4 Installation	14
4.1 Selecting the location	14
4.2 Preparing installation	15
4.3 Installing the instrument	15
5 Operation	16
5.1 Operating controls	16
5.2 Switching on the centrifuge	17
5.3 Replacing the rotor	17
5.3.1 Inserting the rotor	17
5.3.2 Removing the rotor	17
5.4 Closing the centrifuge lid	17
5.5 Preparing for centrifugation	18
5.5.1 Loading the rotor	18
5.5.2 Closing the rotor lid	18
5.5.3 Closing the rotor lid (aerosol-tight centrifugation)	19
5.6 Centrifugation	19
5.6.1 Centrifugation with time setting	19
5.6.2 End of centrifugation	20
5.6.3 Centrifuging in continuous operation	20
5.6.4 Short run centrifugation	20
5.6.5 Setting the alarm volume	20
5.7 Aerosol-tight centrifugation	21
5.8 Switching off the centrifuge	21

Table of contents

6	Maintenance	22
6.1	Service	22
6.2	Prepare cleaning/disinfection	22
6.3	Cleaning/disinfection	23
6.3.1	Cleaning and disinfecting the device	24
6.3.2	Cleaning and disinfecting the rotor	24
6.3.3	Cleaning and disinfecting the rotor lid	24
6.4	Cleaning glass breakage	25
6.5	Replacing fuses	26
6.6	Decontamination before shipment	26
7	Troubleshooting	26
7.1	General errors	26
7.2	Error messages	27
8	Transport, storage and disposal	28
8.1	Transport	28
8.2	Storage	28
8.3	Disposal	28
9	Technical data	29
9.1	Power supply	29
9.2	Ambient conditions	29
9.3	Weight/dimensions	29
9.4	Noise level	29
9.5	Application parameters	30
9.6	Service life of accessories	30
10	Rotors for the Microcentrifuge 24	31
10.1	Rotor Standard 24 AT	31
10.2	Switching on the centrifuge	32
11	Ordering information	33

1 Operating instructions

1.1 Using this manual

- Read this operating manual completely before using the device for the first time. Observe the instructions for use of the accessories where applicable.
- This operating manual is part of the product. Please keep it in a place that is easily accessible.
- Enclose this operating manual when transferring the device to third parties.
- The current version of the operating manual for all available languages can be found on our webpage www.usascientific.com.

1.2 Danger symbols and danger levels

1.2.1 Danger symbols

The safety instructions in this manual have the following danger symbols and danger levels:

	Biohazard		Explosive substances
	Electric shock		Risk of crushing
	Hazard point		Material damage

1.2.2 Danger levels

DANGER	Will lead to severe injuries or death.
WARNING	May lead to severe injuries or death.
CAUTION	May lead to light to moderate injuries.
NOTICE	May lead to material damage.

1.3 Symbols used

Depiction	Meaning
1.	Actions in the specified order
2.	
►	Actions without a specified order
•	List
Text	Display or software texts
	Additional information

1.4 Abbreviations used

PCR = Polymerase Chain Reaction

rcf = Relative centrifugal force: g -force in m/s^2

rpm = Revolutions per minute

UV = Ultraviolet radiation

2 Safety

2.1 Intended use

The Microcentrifuge 24 is used for the separation of aqueous solutions and suspensions of different densities in approved sample tubes.

The Microcentrifuge 24 is exclusively intended for use indoors. All country-specific safety requirements for operating electrical equipment in the laboratory must be observed.

2.2 User profile

The device and accessories may only be operated by trained and skilled personnel. Before using the device, read the operating manual and the instructions for use of the accessories carefully and familiarize yourself with the device's mode of operation.

2.3 Information on product liability

In the following cases, the designated protection of the device may be affected. Liability for any resulting damage or personal injury is then transferred to the owner:

- The device is not used in accordance with the operating manual.
- The device is used outside of its intended use.
- The device is used with accessories or consumables which are not recommended by the manufacturer.
- The device is maintained or repaired by persons who are not authorized by the manufacturer.
- The user makes unauthorized changes to the device.

2.4 Application limits

2.4.1 Declaration concerning the ATEX directive (2014/34/EU)

DANGER! Risk of explosion.

- ▶ Do not operate the device in areas where explosive substances are handled.
- ▶ Do not use this device to process any explosive or highly reactive substances.
- ▶ Do not use this device to process any substances which may generate an explosive atmosphere.

Due to its design and the environmental conditions inside the device, the Microcentrifuge 24 is not suitable for use in a potentially explosive atmosphere. The device may only be used in a safe environment, such as in the open environment of a ventilated laboratory or a fume hood.

The use of substances that may contribute to a potentially explosive atmosphere is not permitted. The final decision on the risks associated with the use of such substances lies with the user.

2.5 Warnings for intended use

2.5.1 Personal injury or damage to device

WARNING! Electric shock due to damage to the device or mains/power cord.

- ▶ Only switch on the device if the device and mains/power cord are undamaged.
- ▶ Only operate devices which have been installed or repaired properly.
- ▶ In case of danger, disconnect the device from the mains/power supply voltage. Disconnect the mains/power plug from the device or the earth/grounded socket. Use the isolating device intended for this purpose (e.g. the emergency switch in the laboratory).

WARNING! Lethal voltages inside the device.

If you touch any parts which are under high voltage you may experience an electric shock. Electric shocks cause injuries to the heart and respiratory paralysis.

- ▶ Ensure that the housing is closed and undamaged.
 - ▶ Do not remove the housing.
 - ▶ Ensure that no liquids can penetrate the device.
- Only authorized service staff may open the device.

WARNING! Danger due to incorrect voltage supply.

- ▶ Only connect the device to voltage sources which correspond with the electrical requirements on the name plate.
- ▶ Only use earth/grounded sockets with a protective earth (PE) conductor.
- ▶ Only use the mains/power cord supplied.

WARNING! Damage to health due to infectious liquids and pathogenic germs.

- ▶ When handling infectious liquids and pathogenic germs, observe the national regulations, the biosafety level of your laboratory, the material safety data sheets, and the manufacturer's application notes.
- ▶ Wear your personal protective equipment.
- ▶ For comprehensive regulations about handling germs or biological material of risk group II or higher, please refer to the "Laboratory Biosafety Manual" (source: World Health Organization, Laboratory Biosafety Manual, the current edition).

WARNING! Risk of injury when opening or closing the centrifuge lid.

There is a risk of crushing your fingers when opening or closing the centrifuge lid.

- ▶ Do not reach between the device and centrifuge lid when opening or closing the centrifuge lid.
- ▶ Do not reach into the locking mechanism of the centrifuge lid.
- ▶ Open the centrifuge lid fully to ensure that the centrifuge lid cannot slam shut.

WARNING! Risk of injury from rotating rotor.

If the emergency release of the lid is operated, the rotor may continue to rotate for several minutes.

- ▶ Wait for the rotor to stop before activating the emergency release.
- ▶ To check, look through the monitoring glass in the centrifuge lid.

WARNING! Risk of injury from chemically or mechanically damaged accessories.

Even minor scratches and cracks can lead to severe internal material damage.

- ▶ Protect all accessory parts from mechanical damage.
- ▶ Inspect the accessories for damage before every use. Replace any damaged accessories.
- ▶ Do not use any accessories which have exceeded their maximum service life.

CAUTION! Poor safety due to incorrect accessories and spare parts.

The use of accessories and spare parts other than those recommended by the manufacturer may impair the safety, functioning and precision of the device. The manufacturer cannot be held liable or accept any liability for damage resulting from the use of accessories and spare parts other than those recommended or from the improper use of such equipment.

- ▶ Only use accessories and original spare parts recommended by the manufacturer.

NOTICE! Damage to the device due to spilled liquids.

1. Switch off the device.
2. Disconnect the device from the mains/power supply.
3. Carefully clean the device and the accessories in accordance with the cleaning and disinfection instructions in the operating manual.
4. If a different cleaning and disinfecting method is to be used, contact the manufacturer to ensure that the intended method will not damage the device.

NOTICE! Damage to electronic components due to condensation.

Condensate may form in the device when it has been transported from a cool environment to a warmer environment.

- ▶ After installing the device, wait for at least 4 h. Only then connect the device to the mains/ power line.

2.5.2 Incorrect handling of the centrifuge

NOTICE! Damage from knocking against or moving the device during operation.

If the rotor hits the rotor chamber wall, it will cause considerable damage to the device and rotor.

- ▶ Do not move or knock against the device during operation.

2.5.3 Incorrect handling of the rotors

WARNING! Risk of injury from improperly attached rotors and rotor lids.

- ▶ Only centrifuge with the rotor and rotor lid firmly tightened.
- ▶ If unusual noises occur when the centrifuge starts, the rotor or rotor lid may not be properly secured. Immediately press the **STOP/OPEN** key to stop centrifuging.

CAUTION! Risk of injury due to asymmetric loading of a rotor.

- ▶ Load rotors symmetrically with identical tubes.
- ▶ Only load adapters with suitable tubes.
- ▶ Always use the same type of tubes (weight, material/density and volume).
- ▶ Check that loading is symmetrical by balancing the adapters and tubes used with a balance.

CAUTION! Risk of injury from overloaded rotor.

The centrifuge is designed for the centrifugation of material with a maximum density of 1.2 g/ml at maximum speed and filling volume and/or load.

- ▶ Do not exceed the maximum load of the rotor.

NOTICE! Damage to rotors from aggressive chemicals.

Rotors are high-quality assemblies which withstand extreme stresses. This stability can be impaired by aggressive chemicals.

- ▶ Avoid using aggressive chemicals such as strong and weak alkalis, strong acids, solutions with mercury ions, copper ions and other heavy metal ions, halogenated hydrocarbons, concentrated saline solutions and phenol.
- ▶ If it is contaminated by aggressive chemicals, clean the rotor and especially the rotor bores immediately using a neutral cleaning agent.

2.5.4 Extreme strain on the centrifugation tubes

CAUTION! Risk of injury from overloaded tubes.

- ▶ Note the loading limits specified by the tube manufacturer.
- ▶ Only use tubes which are approved by the manufacturer for the required *g*-force (rcf).

NOTICE! Risk from damaged tubes.

Damaged tubes must not be used, as this could cause further damage to the device and the accessories and sample loss.

- ▶ Before use, visually check all of the tubes for damage.

NOTICE! Danger due to deformed or brittle tubes. Autoclaving at excessive temperatures can lead to tubes made from plastic becoming brittle and deformed.

This could cause damage to the device and the accessories and sample loss.

- ▶ Observe the temperatures specified by the manufacturer when autoclaving tubes.
- ▶ Do not use brittle or deformed rechargeable tubes.

NOTICE! Danger from open tube lids.

Open tube lids can break off during centrifugation and damage the rotor and the centrifuge.

- ▶ Carefully seal all tube lids before centrifuging.

NOTICE! Damage to plastic tubes from organic solvents.

The strength of plastic tubes is reduced when organic solvents (e.g., phenol, chloroform) are used, which will damage the tubes.

- ▶ Note the manufacturer's information on the chemical resistance of the tubes.

NOTICE! Micro test tubes heat up.

In non-refrigerated centrifuges, the temperature in the rotor chamber, rotor and sample may increase to above 40 °C, depending on the run time, *g*-force (rcf) / speed and ambient temperature.

- ▶ Please note that this will reduce the centrifugation stability of the micro test tubes.
- ▶ Please note the temperature resistance of the samples.

2.6 Safety instructions on device and accessories

Depiction	Meaning	Location
	Caution ▶ Observe the safety instructions in the operating manual.	Right side of the device
	▶ Observe operating manual.	Right side of the device
	Warning of biological risks when handling infectious liquids or pathogenic germs.	Aerosol-tight fixed-angle rotors: Rotor lid

Product description

Microcentrifuge 24

3 Product description

3.1 Product overview


Fig. 3-1: Microcentrifuge 24: Front and side view

- 1 Centrifuge lid**
- 2 Monitoring glass**
For visual control during rotor stop or speed check using a stroboscope
- 3 Control panel**
Display and keys for operating the centrifuge
- 4 Mains/power cord socket**
Connection for the supplied mains/power cord

- 5 Fuse holder**
- 6 Mains/power switch**
Switch for switching the centrifuge on and off
- 7 Name plate**
- 8 Emergency release**

3.2 Delivery package

1	Microcentrifuge 24
1	Rotor key
1	Mains/power cord
1	User manual
1	Set of fuses

- ▶ Check whether the delivery is complete.
- ▶ Check all parts for any transport damage.
- ▶ To safely transport and store the device, retain the transport box and packing material.

3.3 Features

The versatile Microcentrifuge 24 has a capacity of up to 24 × 2 ml and achieves a maximum of 21300 × g or 15060 rpm.

You can select from 2 different rotors to centrifuge the following tubes for various applications:

- Micro test tubes (0.2 ml to 2.0 ml)
- PCR strips
- Microtainers
- Spin columns

3.4 Name plate


Fig. 3-2: Device identification (example)

- | | |
|--|--|
| 1 Maximum density of the material for centrifuging | 8 Maximum rated current |
| 2 Maximum kinetic energy | 9 Maximum rated power |
| 3 Maximum speed | 10 Data matrix code for serial number |
| 4 Serial number | 11 Designation of origin |
| 5 Product name | 12 Certification marks and symbols (device-specific) |
| 6 Rated voltage | 13 Manufacturer and manufacturer's address |
| 7 Rated frequency | |

Tab. 3-1: Certification marks and symbols (device-specific)

Symbol/certification mark	Meaning
	Serial number
	UL listing certification mark: declaration of conformity, USA
	Certification mark for electromagnetic compatibility of the Federal Communications Commission, USA

4 Installation

4.1 Selecting the location

WARNING! Danger due to incorrect voltage supply.

- ▶ Only connect the device to voltage sources which correspond with the electrical requirements on the name plate.
- ▶ Only use earth/grounded sockets with a protective earth (PE) conductor.
- ▶ Only use the mains/power cord supplied.

NOTICE! If an error occurs, any objects in the immediate proximity of the device may become damaged.

- ▶ In accordance with the recommendations of EN 61010-2-020, leave a safety clearance of **30 cm (11.8 in)** around the device during operation.
- ▶ Please remove all materials and objects from this area.

NOTICE! Damage due to overheating.

- ▶ Do not install the device near heat sources (e.g. heating, drying cabinet).
- ▶ Do not expose the device to direct sunlight.
- ▶ Ensure unobstructed air circulation. Maintain a clearance of at least 30 cm (11.8 in) around all ventilation gaps.

NOTICE! Radio interference.

For devices with Class A noise emission in accordance with EN 61326-1/EN 55011, the following applies: This device has been developed and tested in accordance with CISPR 11 Class A. The device may cause radio interference in domestic environments and is not intended for use in residential areas. The device cannot ensure adequate protection of radio reception in residential areas and domestic environments.

- ▶ If necessary, take appropriate measure to eliminate the interferences.

Mains/power connection for centrifuges: The operation of the centrifuge is only permitted in a building installation which complies with the applicable national regulations and standards. In particular, it needs to be ensured that there are no prohibited loads on the supply lines and assemblies that are located before the internal protection of the device. This can be ensured by additional circuit breakers or other suitable fuse elements in the building installation.

The mains/power switch and the disconnecting device of the mains/power line must be easily accessible during operation (e.g. a residual current circuit breaker).

Select the location of the device according to the following criteria:

- Mains/power connection in accordance with the name plate
 - Minimum distance to other devices and walls: 30 cm (11.8 in)
 - Resonance free table with horizontal even work surface
 - The surrounding area must be well ventilated.
 - The location is protected against direct sunlight.
- ▶ Do not use this device near strong electromagnetic sources (e.g. unshielded high frequency sources) as they could impede proper functioning of the device.

4.2 Preparing installation

The weight of the centrifuge is 12.98 kg.

Unpacking the centrifuge

1. Open the packaging board.
2. Remove the accessories.
3. Lift the centrifuge out of the carton.
4. Place the centrifuge on a suitable lab bench.
5. Remove the plastic sleeve.

The device is delivered with and without a rotor.
The transport securing device can be removed directly from devices without a rotor.

6. Turn the rotor nut **counterclockwise** using the supplied rotor key.
7. Lift the rotor out vertically.
8. Remove the transport securing device.

4.3 Installing the instrument

Prerequisites

The device is on a suitable lab bench.

WARNING! Danger due to incorrect voltage supply.

- ▶ Only connect the device to voltage sources which correspond with the electrical requirements on the name plate.
- ▶ Only use earth/grounded sockets with a protective earth (PE) conductor.
- ▶ Only use the mains/power cord supplied.

NOTICE! Damage to electronic components due to condensation.

Condensate may form in the device when it has been transported from a cool environment to a warmer environment.

- ▶ After installing the device, wait for at least 4 h. Only then connect the device to the mains/power line.

1. Let the device warm up to ambient temperature.
2. Connect the centrifuge to the mains and switch it on using the mains/power switch.
 - The **RPM/RCF** key lights up.
 - The display is active.

5 Operation

5.1 Operating controls


Fig. 5-1: Microcentrifuge 24 operating controls

1 Display

2 START key

Start centrifugation

3 STOP/OPEN key

Stop the centrifugation and release the lid

4 SPEED arrow keys

Set the speed of centrifugation

Touch and hold the arrow key: quick setting

5 RPM/RCF key

Switch the display of the centrifugation speed (rpm or rcf)

6 SHORT key

Short run centrifugation

7 TIME arrow keys

Set the centrifugation time

Touch and hold the arrow key: quick setting


Fig. 5-2: Display Microcentrifuge 24

1 Centrifugation time

2 g-force (rcf) or rotational speed (rpm)

Actual value

5.2 Switching on the centrifuge

- ▶ Switch the centrifuge on at the mains/power switch.
 - The parameter settings of the last run are displayed.
 - The lid opens.

5.3 Replacing the rotor

5.3.1 Inserting the rotor

1. Place the rotor vertically onto the motor shaft from the top.
2. Insert the supplied rotor key into the rotor nut.
3. Turn the rotor key **clockwise** until the rotor nut is firmly tightened.

5.3.2 Removing the rotor

1. Turn the rotor nut **counterclockwise** using the supplied rotor key.
2. Remove the rotor by lifting it vertically.

5.4 Closing the centrifuge lid

- ▶ Check correct attachment of rotor and rotor lid.
- ▶ Push down the centrifuge lid until the lid latch engages and the lid is automatically closed.

The centrifuge will close automatically.

The **STOP/OPEN** key lights up blue.

5.5 Preparing for centrifugation

5.5.1 Loading the rotor

CAUTION! Risk of injury due to asymmetric loading of a rotor.

- ▶ Load rotors symmetrically with identical tubes.
- ▶ Only load adapters with suitable tubes.
- ▶ Always use the same type of tubes (weight, material/density and volume).
- ▶ Check that loading is symmetrical by balancing the adapters and tubes used with a balance.

1. Check the maximum payload (adapter, tube and contents) for each rotor bore.
2. Load rotors and adapters only with the tubes intended for them.
3. To ensure symmetrical loading, insert sets of two tubes in opposite bores. Tubes located opposite each other must be of the same type and contain the same filling quantity.


Fig. 5-3: Symmetrical loading of a fixed-angle rotor

To keep the weight differences between the filled tubes low, we recommend taring with a balance. This is gentle on the drive and reduces the operating noises.

5.5.2 Closing the rotor lid

Use matching rotor lids

- Fixed-angle rotors may only be used with the appropriate rotor lid for the respective rotor. The rotor name on the rotor must correspond to the rotor name on the rotor lid.

1. Place the rotor lid vertically on the rotor.
2. Turn the rotor lid screw clockwise to seal the rotor.

5.5.3 Closing the rotor lid (aerosol-tight centrifugation)

i Identification of aerosol-tight rotors


An aerosol-tight rotor and the matching aerosol-tight rotor lid must be used for aerosol-tight centrifugation.

Aerosol-tight fixed-angle rotor

- Designation ends with **AT**

Aerosol-tight rotor lid

- Labeled **aerosol-tight**

1. Check the correct positioning of the external sealing ring in the groove.
2. Place the rotor lid on the rotor in a vertical motion.
3. Turn the rotor lid screw clockwise until it stops to seal the rotor.

5.6 Centrifugation

Prerequisites

- The centrifuge is switched on.
- The rotor has been inserted and attached correctly.
- The rotor has been loaded correctly.
- The rotor lid has been mounted correctly.
- The centrifuge lid is closed.

WARNING! Risk of injury from improperly attached rotors and rotor lids.

- ▶ Only centrifuge with the rotor and rotor lid firmly tightened.
- ▶ If unusual noises occur when the centrifuge starts, the rotor or rotor lid may not be properly secured. Immediately press the **STOP/OPEN** key to stop centrifuging.

5.6.1 Centrifugation with time setting

Setting the centrifugation parameters

1. Set the centrifugation time with the **TIME** arrow keys.
2. Set the rotational speed (rpm) or *g-force* (rcf) with the **SPEED** arrow keys.

Starting the centrifugation run

3. To start the centrifugation run, press the **START** key.

Switch the centrifugation speed display using the **RPM/RCF** key.

Display during centrifugation

- The  symbol is flashing when the rotor is running.
- Remaining run time in minutes. The last minute is counted down in seconds.
- Current *g-force* (rcf) and/or rotational speed (rpm).

During the run you can change the following parameters:

- Centrifugation time
- Speed: During the run, you can switch between the *g-force* and the rotational speed (rpm) using the **RPM/RCF** key.

The following keys are blocked during centrifugation:

- **SHORT** key

5.6.2 End of centrifugation

- ▶ Press the **STOP/OPEN** key to end centrifugation before the set time.
 - After completion of the set time, the centrifuge stops automatically.
 - After the run, the lid opens automatically.
 - During the braking process, the elapsed running time flashes on the display.
 - If the speaker is switched on, a signal sounds when the rotor has stopped.

5.6.3 Centrifuging in continuous operation

Setting continuous run

1. In order to centrifuge without any time limits, use the **TIME** arrow keys to select the setting ∞ (▼ below 10 s or ▲ above 9:59 h).
2. Set the rotational speed (rpm) or *g-force* (rcf) with the **SPEED** arrow keys.
3. To start the centrifugation run, press the **START** key.
 - The  symbol is flashing when the rotor is running.
 - The cycle time is counted up.
 - Current *g-force* (rcf) and/or rotational speed (rpm).

5.6.4 Short run centrifugation

Short run centrifugation takes places up to the maximum rotational speed of the rotor used (15060 rpm).

- ▶ Keep the **SHORT** key pressed.
The centrifuge accelerates up to the maximum rotational speed of the rotor used and continues to run at maximum speed.
- ▶ Release the **SHORT** key.
The centrifuge stops.

5.6.5 Setting the alarm volume

Prerequisites

The lid is open.

- ▶ Touch and hold the **START** key and set the desired volume using the **TIME** arrow keys.

5.7 Aerosol-tight centrifugation

**WARNING! Risk to health due to limited aerosol tightness with incorrect rotor/rotor lid combination.**

Aerosol-tight centrifugation is guaranteed only if the rotors and rotor lids intended for this purpose are used. The designation of an aerosol-tight fixed-angle rotor ends with **AT**.

- ▶ Always use a rotor and rotor lid marked aerosol-tight together for aerosol-tight centrifugation. The details specifying in which centrifuge you may use the aerosol-tight rotor and rotor lid can be found on the rotor and on the top of the rotor lid.
- ▶ Only use the aerosol-tight rotor lid in combination with the rotor that is specified on the rotor lid.

**WARNING! Damage to health as a result of limited aerosol tightness and incorrect usage.**

Mechanical stresses and contamination by chemicals or other aggressive solvents may impair the aerosol tightness of the rotors and rotor lids. Autoclaving at excessive temperatures can lead to vessels, adapters and rotor lids becoming brittle and deformed.

- ▶ Check the integrity of the seals of the aerosol-tight rotor lids or caps before each use.
- ▶ Only use aerosol-tight rotor lids or caps if the seals are undamaged and clean.
- ▶ Do not exceed temperatures of 121°C or a time of more than 20 min. while autoclaving.
- ▶ Lightly grease the threads of the rotor lid screw with pivot grease after every proper autoclaving process (121 °C, 20 min.).
- ▶ For aerosol-tight rotor lids with exchangeable seal, the seal must be replaced after 50 autoclaving cycles.
- ▶ **Never** store aerosol-tight rotors or buckets closed.

The aerosol tightness of rotors, rotor lids, buckets and caps has been tested and certified according to Annex AA of IEC 61010-2-020

5.8 Switching off the centrifuge

1. Open the centrifuge lid.
Residual moisture can evaporate.
2. Remove rotor lids from fixed-angle rotors.
Aerosol-tight accessories may not be stored with the lid closed.
3. Switch off the centrifuge using the mains/power switch.

6 Maintenance
6.1 Service

WARNING! Risk of fire or electrical shock

- ▶ Have the centrifuge's electrical safety, especially the paths for the protective connections, checked every 12 months by trained and skilled personnel.

Please note the country-specific regulations.

6.2 Prepare cleaning/disinfection

- ▶ Clean all accessible surfaces of the device and the accessories at least weekly and when contaminated.
- ▶ Clean the rotor regularly. This way the rotor is protected and the durability is prolonged.
- ▶ Furthermore, observe the notes on decontamination (see *Decontamination before shipment on p. 26*) when the device is sent to the authorized Technical Service for repairs.

The procedure described in the following chapter applies to the cleaning as well as to the disinfection or decontamination. The table below describes the steps required on top of this:

Cleaning	Disinfecting/decontamination
1. Use a mild cleaning fluid to clean the accessible surfaces of the device and the accessories. 2. Carry out the cleaning as described in the following chapter.	1. Choose the disinfection method which corresponds to the legal regulations and guidelines in place for your range of application. For example, use alcohol (ethanol, isopropanol) or alcohol-based disinfectants. 2. Carry out the disinfection or decontamination as described in the following chapter. 3. Then clean the device and the accessories.

6.3 Cleaning/disinfection

DANGER! Electric shock due to the ingress of liquid.

- ▶ Switch off the device and disconnect it from the mains/power line before starting cleaning or disinfection.
- ▶ Do not allow any liquids to penetrate the inside of the housing.
- ▶ Do not perform a spray clean/spray disinfection on the housing.
- ▶ Only reconnect the device to the mains/power line when it is completely dry, both inside and outside.

WARNING! Damage to health as a result of limited aerosol tightness and incorrect usage.

Mechanical stresses and contamination by chemicals or other aggressive solvents may impair the aerosol tightness of the rotor and rotor lid. Autoclaving at excessive temperatures can lead to vessels, adapters and rotor lids becoming brittle and deformed.

- ▶ Check the integrity of the seal of the aerosol-tight rotor lid before each use.
- ▶ Only use aerosol-tight rotor lids with undamaged and clean seals.
- ▶ Do not exceed temperatures of 121°C or a time of more than 20 min. while autoclaving.
- ▶ Lightly grease the threads of the rotor lid screw with pivot grease after every proper autoclaving process (121 °C, 20 min.).
- ▶ **Never** store the aerosol-tight rotor or bucket closed.

NOTICE! Damage from the use of aggressive chemicals.

- ▶ Do not use any aggressive chemicals on the device or its accessories, such as strong and weak bases, strong acids, acetone, formaldehyde, halogenated hydrocarbons or phenol.
- ▶ If the device has been contaminated by aggressive chemicals, clean it immediately using a mild cleaning agent.

NOTICE! Corrosion due to aggressive cleaning agents and disinfectants.

- ▶ Do not use any corrosive cleaning agents, aggressive solvents or abrasive polishes.
- ▶ Do not incubate the accessories in aggressive cleaning agents or disinfectants for longer periods.

NOTICE! Damage from UV and other high-energy radiation.

- ▶ Do not use UV, beta, gamma, or any other high-energy radiation for disinfection.
- ▶ Avoid storage in areas with strong UV radiation.

i Autoclaving

Rotors, rotor lids and adapters can be autoclaved (121 °C, 20 min).
Replace the seal of the aerosol-tight rotor lid after 50 autoclaving cycles.

i Aerosol tightness

Check that the seal is intact before use.
Replace the rotor lid with screw cap when the sealing ring on the lid screw and in the lid groove become worn.
Regular care of the sealing ring is necessary in order to protect the rotor.
Never store the aerosol-tight rotor with the lid screwed on!
In order to prevent damage, lightly grease the lid thread of the aerosol-tight rotor with pivot grease.

6.3.1 Cleaning and disinfecting the device

1. Open the lid. Switch the device off at the mains/power switch. Disconnect the mains/power plug from the voltage supply.
2. Remove the rotor.
3. Clean and disinfect all accessible surfaces on the device including the mains/power cord using a damp cloth and recommended cleaning agents.
4. Thoroughly clean the rubber seal of the rotor chamber with water.
5. Rub the dry rubber seal with glycerol or talcum powder to prevent it from becoming brittle. Other components of the device, such as the motor shaft and rotor cone, must not be lubricated.
6. Clean the motor shaft with a soft, dry, lint-free cloth. Do not grease the motor shaft.
7. Check the motor shaft for damage.
8. Check the device for corrosion and damage.
9. Leave the centrifuge lid open when the device is not being used.
10. Only connect the device to the power supply if it is fully dry inside and out.

6.3.2 Cleaning and disinfecting the rotor

1. Inspect the rotor and accessories for damage and corrosion. Do not use any damaged rotors or accessories.
2. Clean and disinfect the rotors and accessories using the recommended cleaning agents.
3. Clean and disinfect the rotor bores using a bottle brush.
4. Rinse the rotors and accessories thoroughly with distilled water. Rinse the rotor bores of fixed-angle rotors particularly thoroughly.

i Do not immerse the rotor in liquid as liquid can enter through the openings when doing so.

5. Place the rotors on a towel to dry. Place fixed-angle rotors with the rotor bores facing downwards to allow the bores to also dry.
6. Clean the rotor cone with a soft, dry, lint-free cloth. Do not lubricate the rotor cone.
7. Inspect the rotor cone for damage.
8. Place the dry rotor onto the motor shaft.
9. Tighten the rotor nut by turning it **clockwise**.
10. Leave the rotor lid open when the rotor is not being used.


6.3.3 Cleaning and disinfecting the rotor lid

Prerequisites


The rotor lid has been removed.

Recommended cleaning agents:

- Alcohol 70% (ethanol, isopropanol)
- Mild, neutral cleaning agent

1. Remove the sealing ring to thoroughly clean the groove below it.
2. Clean and disinfect the rotor lid using the recommended cleaning agents.
3. Rinse the rotor lid thoroughly with distilled water.

4. Moisten the new sealing ring with clean water.
5. Insert the sealing ring in the clean groove of the rotor lid.
6. Press the sealing ring into the lateral groove, around the entire circumference of the rotor lid.
7. Place the rotor lid with the underside facing upwards on a cloth.
8. Leave the rotor lid to dry for 5 –10 minutes.
9. Perform a visual inspection.
The seal must be flush with the groove of the rotor lid around the entire circumference and must not protrude at any point.
10. Fit the rotor lid on the rotor.
11. Leave the rotor lid open when the rotor is not being used.

i The rotor lid cannot close properly if the sealing ring is not correctly inserted.

6.4 Cleaning glass breakage

When using glass tubes there is a risk of glass breakage in the rotor chamber. The resulting glass splinters are swirled around in the rotor chamber during centrifugation and have a sandblasting effect on the rotor and accessories. Smallest glass particles become lodged in the rubber parts (e.g., the motor sleeve, the rotor chamber seal, and the rubber mats of adapters).

NOTICE! Glass breakage in the rotor chamber

Glass tubes in the rotor chamber may break if the *g*-force is too high. Broken glass can damage the rotor, accessories and samples.

- Please note the manufacturer's information on the recommended centrifugation parameters (load and speed).

Effects of glass breakage in the rotor chamber:

- Fine black metal abrasion dust in the rotor chamber (in metal rotor bowls).
- The surfaces of the rotor chamber and accessories are scratched.
- The chemical resistance of the rotor chamber is reduced.
- Contamination of samples.
- Wear on rubber parts.

How to proceed in case of glass breakage

1. Remove all splinters and glass powder from the rotor chamber and accessories.
2. Thoroughly clean the rotor and rotor chamber. Thoroughly clean the bores of the fixed-angle rotors, in particular.
3. Regularly check the rotor bores for deposits and damage.

6.5 Replacing fuses

The fuse holder is located under the mains power socket.

1. Disconnect the mains/power plug.
2. Open the fuse holder.
3. Replace faulty fuses and reinsert the fuse holder.

6.6 Decontamination before shipment

Please contact us before shipping a device.

If you are shipping the device to the authorized Technical Service for repairs or to your authorized dealer for disposal please note the following:

WARNING! Risk to health from contaminated device.

1. Observe the information on the decontamination certificate. It is available as a PDF document on our webpage www.usascientific.com.
2. Decontaminate all the parts you are going to dispatch.
3. Include the fully completed decontamination certificate in the shipment.

7 Troubleshooting

If you are unable to resolve the error with the suggested measures, please contact your local sales representative.

7.1 General errors

Problem	Cause	Solution
No display.	No mains/power connection. Mains/power outage.	▶ Check the mains/power connection. ▶ Check the fuse of the device. ▶ Check the mains/power fuse of the lab.
Centrifuge lid cannot be opened.	The rotor is still running.	▶ Wait for the rotor to stop.
	Mains/power outage.	1. Check the fuse of the device. 2. Check the mains/power fuse of the lab. 3. Actuate the emergency lid release.
Centrifuge cannot be started.	The centrifuge lid is not closed.	▶ Close the centrifuge lid.
Centrifuge shakes when it starts up.	The rotor is loaded asymmetrically.	1. Stop the centrifuge and load the rotor symmetrically. 2. Restart the centrifuge.

7.2 Error messages

If an error message appears, proceed as follows:

1. Remedy the fault as described in the "Solution" column.
2. To clear the error message from the display, press the **STOP/OPEN** key.
3. If necessary, repeat centrifugation.

Problem	Cause	Solution
IMBAL	The rotor is loaded asymmetrically.	► Load the rotor symmetrically and balance it.
NET INT	Mains/power failure during a run.	► Check the mains/power supply.
LID ERROR	Centrifuge lid cannot be locked.	► Try to close the centrifuge lid again.
	Centrifuge lid cannot be released.	1. Switch off the centrifuge and wait for 20 s. 2. Switch on the centrifuge. If the error occurs again: 1. Switch off the centrifuge. 2. Actuate the emergency lid release.
	Prohibited opening of lid during a run or lid switch defective.	1. Wait for the rotor to stop. 2. Open the centrifuge lid and then close it again. 3. Repeat the run.
LID LIFT	The centrifuge lid has not been opened wide enough.	► Open the centrifuge lid wider by hand.
NO RPM	Error in the rotational speed measurement system.	► Leave the device switched on until the rotor stops and the error message disappears (up to 15 min).
ERROR 6	Error in the drive electronics.	► Repeat the run. If the error message appears again: 1. Switch off the centrifuge and wait for 20 s. 2. Switch on the centrifuge.
ERROR 7	Deviation in the speed check.	1. Wait for the rotor to stop. 2. Tighten the rotor.
ERROR 10	Error during initialization or in the memory.	1. Switch off the centrifuge and wait for 20 s. 2. Switch on the centrifuge.
ERROR 16	Data communication error with the operation control.	1. Switch off the centrifuge and wait for 20 s. 2. Switch on the centrifuge.
ERROR 20	Drive overheated.	► Allow the drive to cool down for at least 15 min.
ERROR 26	Data communication error with the motor.	1. Switch off the centrifuge and wait for 20 s. 2. Switch on the centrifuge.
ERROR 27	Electronics fault.	1. Switch off the centrifuge and wait for 20 s. 2. Switch on the centrifuge.

8 Transport, storage and disposal

8.1 Transport

- ▶ Remove the rotor from the centrifuge before transport.
- ▶ Use the original packing for transport.

	Air temperature	Relative humidity	Atmospheric pressure
General transport	- 25 °C – 60 °C	10 % – 75 %	30 kPa – 106 kPa
Air freight	- 20 °C – 55 °C	10 % – 75%	30 kPa – 106 kPa

8.2 Storage

	Air temperature	Relative humidity	Atmospheric pressure
In transport packing	- 25 °C – 55 °C	10 % – 75 %	70 kPa – 106 kPa
Without transport packing	- 5 °C – 45 °C	10 % – 75 %	70 kPa – 106 kPa

8.3 Disposal

If the product needs to be disposed of, the relevant legal regulations must be observed.

9 Technical data

9.1 Power supply

Microcentrifuge 24

Mains/power connection	120 V, 50 Hz – 60 Hz
Current consumption	4.0 A
Power consumption	270 W
EMC: noise emission (radio interference)	CFR 47, Part 15 – Class B
EMC: noise immunity	EN 61326-1
Overvoltage category	II
Protection class	I / 250V 8AT HBC
Fuses – 120 V	8.0 A
Degree of pollution	2

9.2 Ambient conditions

Environment	For indoor use only
Ambient temperature	2 °C – 40 °C
Relative humidity	10 % – 80 %, non-condensing
Atmospheric pressure	75 kPa – 106 kPa

9.3 Weight/dimensions

Dimensions	Width: 24.2 cm / 9.5 in Depth: 34.6 cm / 13.6 in Height: 21.9 cm / 8.6 in
Weight without rotor	12.98 kg
Rotor weights:	
Standard 24 AT	772 g
PCR 4 × 8	460 g

9.4 Noise level

The noise level was measured according to DIN EN ISO 3745 frontally in a sound measuring room with accuracy class 1 at a distance of 1 m from the device and at lab bench height.

Noise level	< 56 dB(A)
-------------	------------

9.5 Application parameters

Tab. 9-1: Acceleration time and braking time according to DIN 58 970

Rotor	Acceleration time	Deceleration time
Standard 24 AT	15 s	15 s
PCR 4x8	15 s	15 s
Cycle time	10 s – 9:59 h, unlimited (∞) • 10 s – 2 min: can be set in increments of 10 s • 2 min – 10 min: can be set in increments of 30 s • 10 min – 9:59 h: can be set in increments of 1 min	
Rotational speed	100 rpm – 15060 rpm • 100 rpm – 5000 rpm: can be set in increments of 10 rpm • 5000 rpm – 15060 rpm: can be set in increments of 100 rpm	
Relative centrifugal force	10 $\times g$ – 21300 $\times g$ • 10 $\times g$ – 3000 $\times g$: can be set in increments of 10 $\times g$ • 5000 $\times g$ – 21300 $\times g$: can be set in increments of 100 $\times g$	
Maximum load	Fixed-angle rotor: 24 $\times$ 2 ml	
Maximum kinetic energy	4.12 kJ	
Permitted density of the material for centrifuging (at maximum <i>g-force</i> (rcf) or rotational speed (rpm) and maximum load)	1.2 g/ml	

9.6 Service life of accessories

CAUTION! Danger due to material fatigue.

When the service life is exceeded, it cannot be guaranteed that the material of the rotors and accessories will withstand the stresses during centrifugation.

- Do not use any accessories which have exceeded their maximum service life.

All rotors and rotor lids can be used during the entire service life of the centrifuge if the following conditions are met:

- proper use
- recommended maintenance
- undamaged condition

Accessories	Maximum service life after initial setup	
Aerosol-tight rotor lid	–	3 years
Seals in the aerosol-tight rotor lid	50 autoclaving cycles	–
Adapter	–	1 year

The date of manufacture is stamped on the rotors in the format 2015-03 (= March 2015).


10 Rotors for the Microcentrifuge 24

Only use rotors that are intended for use with the corresponding centrifuge.

Please note the manufacturer's information on the centrifugation resistance of the sample tubes used (maximum *g-force*).

10.1 Rotor Standard 24 AT

Aerosol-tight fixed-angle rotor for 24 tubes

	Max. <i>g-force</i> :	21300 × <i>g</i>
	Max. rotational speed:	15060 rpm
Rotor Standard 24 AT	Max. load (adapter, tube and contents):	24 × 3.75 g

Tube	Tube Capacity Tubes per adapter/rotor	Adapter Order no. (international)	Bottom shape Diameter	Max. <i>g-force</i> Max. rotational speed Radius
	PCR tube 0.2 ml 1/24	 8030-0014	Conical Ø 6 mm	21300 × <i>g</i> 15060 rpm 6.3 cm
	Micro test tube 0.4 ml 1/24	 8030-0015	Conical Ø 6 mm	21300 × <i>g</i> 15060 rpm 8.4 cm
	Micro test tube 0.5 ml 1/24	 8030-0013	– Ø 8 mm	21300 × <i>g</i> 15060 rpm 7.3 cm
	Microtainers 0.6 ml 1/24	 8030-0013	– Ø 8 mm	21300 × <i>g</i> 15060 rpm 8.4 cm
	Micro test tube 1.5 ml/2 ml –/24	–	Conical Ø 11 mm	21300 × <i>g</i> 15060 rpm 8.4 cm


Rotors for the Microcentrifuge 24

Microcentrifuge 24

10.2 Rotor PCR 4×8

Fixed-angle rotor for PCR strips and PCR tubes

	Max. <i>g</i> -force:	18257 × <i>g</i>
	Max. rotational speed:	15060 rpm
Rotor PCR 4×8	Max. load (adapter, tube and contents):	4 × 3.5 g

Tube	Tube Capacity Vessels per rotor	Bottom shape Diameter	Max. <i>g</i> -force Max. rotational speed Radius
	PCR strips 8 × 0.2 ml 4 × 8	Conical Ø 6 mm	18257 × <i>g</i> 15060 rpm 7.2 cm
	PCR tube 0.2 ml 32	Conical Ø 6 mm	18257 × <i>g</i> 15060 rpm 7.2 cm

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