


Easy-to-Remove Polycarbonate Lid



Press on the tabs to release the lid

Micronic polycarbonate ComoRack-96 and ComoRack-24 racks are now available with easy-to-remove polycarbonate lids. Just press on the tabs at each end to quickly release the lid locks. The dimensions of the new lid are the same as the old lid, and lids can be purchased separately for use with the original Micronic polycarbonate racks.

Catalog No.	Description	Quantity
1776-2012	New, easy-to-remove polycarbonate lid	10 lids
1772-2502	1.4 ml U-bottom non-coded tubes in ComoRack-96	10 racks of 96
1772-2501	1.4 ml U-bottom alphanumeric tubes in ComoRack-96	10 racks of 96
1775-2311	1.4 ml U-bottom Tracker 2D tubes in ComoRack-96	10 racks of 96
1772-2504	Polycarbonate ComoRack-96 for 1.4 ml tubes (w/o tubes)	10 racks
1772-3100	Polycarbonate ComoRack-24 for 7.5 ml tubes (w/o tubes)	10 racks
1773-2029	7.5 ml non-coded tubes in polycarbonate ComoRack-24	10 racks of 24



MICRONIC

Traceable Sampling Solutions

Please see www.usascientific.com for more information on Micronic tubes, capping options, scanners, and tube sorters.