


# Micronic Manual Capmatsealer

**Designed to ensure correct and uniform application of Micronic Capmats and Capclusters on tubes in racks.**



The Micronic Manual Capmatsealer with an adapter plate for racks of 0.50 ml tubes

## Benefits of Micronic Manual Capmatsealers:

- Compatible with Micronic tubes, racks, TPE Caps/Capclusters/Capbands/Capmats, and EVA Capmats/Capbands
- Caps 1.4 ml, 0.50 ml, 0.75 ml, and 1.1 ml Micronic tubes (*1.4 ml is standard - smaller tubes require adapter blocks*)
- Evenly and correctly caps tubes in racks
- Caps entire rack with one slide of the sealing handle
- Body of Capmatsealer constructed of stainless steel and anodized aluminum
- Body (excluding feet) is autoclavable

Item #	Description
1765-5006	Manual Capmatsealer for 1.4 ml, 0.50 ml, and 0.75 ml tubes ( <i>1.4 ml is standard; adapters are required for smaller tubes</i> )
1765-5702	Support for 0.50 ml tube rack
1765-5700	Support for 0.75 ml tube rack
1765-5701	Support for 1.1 ml tube rack



Please see [www.usascientific.com](http://www.usascientific.com) for more information on Micronic tubes, capping options, scanners, and tube sorters.