


NEW! Automated Screw Cap Handler



Decap or recap an entire rack of 96 screw cap tubes in less than 15 seconds. User-friendly unit saves time and reduces the chance of repetitive strain injury.

- Processes all 96 tubes simultaneously
- Ensures each cap receives the same amount of torque for optimal seal quality
- Can be used in a laboratory or environmental control cabinet
- Does not require a PC interface or control
- Anti-collision design ensures a new rack can not be loaded when the unit is already holding caps



Bench top unit is made from solvent resistant materials. Measures 14.3" W x 11.7" D x 16.1" H. Weighs under 100 pounds.

1778-5200



Please see www.usascientific.com for more information on Micronic tubes, capping options, scanners, and tube sorters.