

# Micronic Manual Decappers

Designed for quick and easy removal of TPE Caps



The Decapper-1 removes single TPE caps from racks of tubes quickly and easily



The Aluminum Decapper-8 was created after a customer requested a decapper sturdy enough to handle frozen tubes




The Decapper-8 removes a row of 8 TPE caps in one motion

## Benefits of Micronic Manual Decappers:

- Designed to remove TPE caps individually or in rows of eight at a time
- Quick and convenient removal - caps are removed smoothly each time
- Aluminum Decapper-8 is sturdy enough to withstand decapping frozen tubes
- Clean room (Class 10,000) production

Item #	Description
1765-4000	Manual Decapper-1, each
1765-4001	Manual Decapper-8, each
1765-5010	Manual Decapper-8 Aluminum, each



**MICRONIC**

Traceable Sampling Solutions

Please see [www.usascientific.com](http://www.usascientific.com) for more information on Micronic tubes, capping options, Scanners, and tube sorters