


NEW Micronic Electric Capper

A fast, easy-to-use, semi-automated capmatsealer for capping a variety of tube sizes with capmats and capclusters.



The New Electric Capper quickly and evenly caps racks of tubes arranged in a 96-well format.

Benefits of the New Electric Capper:

- Versatile: caps all Micronic tubes using capmats and capclusters
*Also compatible with other brands of tubes
- Easy-to-use:
 - *Place rack of tubes inside drawer
 - *Place caps on rack
 - *Close drawer
 - *Press the button to cap tubes
- Flexible: caps both full and partially-full racks of tubes
- The most secure seal: controlled-pressure system caps each tube evenly
- Fast: caps entire rack of 96 tubes in 5 seconds

Features of the Electric Capper:

- Compatible with the Micronic capping range, tubes and racks
- Also compatible with other brands of tubes and caps
(Please contact us for compatible brands)
- Robust, durable bench-top device
- 1.40 ml tubes standard - adapters necessary for smaller tubes
- Power: 230/115V-1
- Weight and dimensions: 29 lbs, 15" x 6" x 13"



Item #	Description
1763-5015	Micronic Electric Capper
1765-5702	Adapter for 0.50 ml Screw Cap tubes
1765-5700	Adapter for 0.75 ml tubes
1765-5701	Adapter for 1.10 ml Screw Cap tubes



Please see www.usascientific.com for more information on Micronic tubes, capping options, scanners, and tube sorters.

Micronic Capping Range



2008 Cat. #	TPE (Thermo Plastic Elastomer) Caps	Pack Size
1775-3003	Pierceable TPE Cap for capping 1 tube bulk	960 pieces
	Capmats	
1775-3000	Pierceable TPE Capmat-96 Natural	10 pieces
1775-3011	Pierceable TPE Capmat-96 Blue	10 pieces
1775-3012	Pierceable TPE Capmat-96 Green	10 pieces
1775-3013	Pierceable TPE Capmat-96 Grey	10 pieces
1775-3014	Pierceable TPE Capmat-96 Pink	10 pieces
1775-3015	Pierceable TPE Capmat-96 Brown	10 pieces
1775-3016	Pierceable TPE Capmat-96 Black	10 pieces
1775-3017	Pierceable TPE Capmat-96 Yellow	10 pieces
	Capbands	
1775-3002	Pierceable TPE Capband-8 Natural in Capmat form	10 pieces
1775-3018	Pierceable TPE Capband-8 Blue in Capmat form	10 pieces
1775-3019	Pierceable TPE Capband-8 Green in Capmat form	10 pieces
1775-3020	Pierceable TPE Capband-8 Grey in Capmat form	10 pieces
1775-3021	Pierceable TPE Capband-8 Pink in Capmat form	10 pieces
1775-3022	Pierceable TPE Capband-8 Brown in Capmat form	10 pieces
1775-3023	Pierceable TPE Capband-8 Black in Capmat form	10 pieces
1775-3024	Pierceable TPE Capband-8 Yellow in Capmat form	10 pieces
1775-3035	Pierceable TPE Capband-12 Natural in Capmat form <i>(Please enquire about additional colors)</i>	10 pieces
	Capclusters	
1775-3001	Pierceable TPE Capcluster-96 Natural	10 pieces
1775-3004	Pierceable TPE Capcluster-96 Blue	10 pieces
1775-3005	Pierceable TPE Capcluster-96 Green	10 pieces
1775-3006	Pierceable TPE Capcluster-96 Grey	10 pieces
1775-3007	Pierceable TPE Capcluster-96 Pink	10 pieces
1775-3008	Pierceable TPE Capcluster-96 Rust	10 pieces
1775-3009	Pierceable TPE Capcluster-96 Black	10 pieces
1775-3010	Pierceable TPE Capcluster-96 Yellow	10 pieces
1775-3036	Pierceable TPE Capcluster-96 White	10 pieces
1775-3038	Pierceable TPE Capcluster-96 Purple	10 pieces
1775-3039	Pierceable TPE Capcluster-96 Light Blue	10 pieces
1775-3037	Pierceable TPE Capcluster-96 in 8 colors: (Natural, Blue, Green, Grey, Pink, Brown, Black and Yellow)	10 pieces
	Split Caps	
1775-3027	Split TPE Capband-8	10 pieces
1775-3026	Split TPE Capcluster-96 Natural	10 pieces
1775-3028	Split TPE Capcluster-96 Blue	10 pieces
1775-3029	Split TPE Capcluster-96 Green	10 pieces
1775-3030	Split TPE Capcluster-96 Grey	10 pieces
1775-3031	Split TPE Capcluster-96 Pink	10 pieces
1775-3032	Split TPE Capcluster-96 Brown	10 pieces
1775-3033	Split TPE Capcluster-96 Black	10 pieces
1775-3034	Split TPE Capcluster-96 Yellow	10 pieces
	Accessories	
1765-4000	Manual Decapper-1	1 piece
1765-4001	Manual Decapper-8	1 piece
1765-5010	Manual Decapper-8 Aluminum	1 unit
1763-5000	Automatic Decapper-96	1 unit