

NEW! Micronic TPE Cap Range... Seals better than ever!

Micronic offers the widest range of TPE (Thermo Plastic Elastomer) caps - 11 colors are available. TPE caps have also been improved with a new elastomer formulation.



Enhanced Elastomer Benefits:

- Better sealing
- More hydrophobic with low permeability to moisture, air and oxygen
- Easier to pierce; reseals better after piercing
- Medical approved material (USP Class VI certificate)
- Excellent compatibility with other Micronic products

TPE Cap Benefits:

- May be pierced multiple times without losing the seal integrity
- Working temperature -80C to +80C
- Split septum option can be pierced with a pipet tip
- Produced in a Class 10,000 clean room
- RNase and DNase free
- Provides superior chemical resistance compared to standard EVA and silicon capmats
- Available sterile (special order)

TPE Cap Range includes:


- **Single caps for capping individual tubes.**
- **Capband-8 and -12 in Capmat format.** Fixed-cap mat holds 96 TPE caps in perforated strips of 8 or 12 for easy separation into strips.
- **Capmat.** Fixed-cap mat seals the entire rack of 96 tubes. It can be removed in one piece and reused.
- **Capcluster.** Seals 96 tubes and has a supporting film that is removed once the rack is sealed to create individually capped tubes.
- **Split cap.** Reduces the piercing force required. May be used with needles, probes, or pipet tips.



MICRONIC

Traceable Sampling Solutions

Cap Removal Tools:



Manual Decapper-8



Manual Capmatsealer



Pneumatic Capper



Automatic Decapper

CATALOG NO.	TPE (Thermo Plastic Elastomer) Caps	QUANTITY
1775-3003	Piercable TPE Cap for capping 1 tube, bulk	960 pieces
	Capmats	
1775-3000	Piercable TPE Capmat Natural for capping 96 tubes	10 mats
1775-3011	Piercable TPE Capmat Blue for capping 96 tubes	10 mats
1775-3012	Piercable TPE Capmat Green for capping 96 tubes	10 mats
1775-3013	Piercable TPE Capmat Grey for capping 96 tubes	10 mats
1775-3014	Piercable TPE Capmat Pink for capping 96 tubes	10 mats
1775-3015	Piercable TPE Capmat Brown for capping 96 tubes	10 mats
1775-3016	Piercable TPE Capmat Black for capping 96 tubes	10 mats
1775-3017	Piercable TPE Capmat Yellow for capping 96 tubes	10 mats
	Capbands	
1775-3002	Piercable TPE Capband-8 Natural in Capmat format	10 mats
1775-3018	Piercable TPE Capband-8 Blue in Capmat format	10 mats
1775-3019	Piercable TPE Capband-8 Green in Capmat format	10 mats
1775-3020	Piercable TPE Capband-8 Grey in Capmat format	10 mats
1775-3021	Piercable TPE Capband-8 Pink in Capmat format	10 mats
1775-3022	Piercable TPE Capband-8 Brown in Capmat format	10 mats
1775-3023	Piercable TPE Capband-8 Black in Capmat format	10 mats
1775-3024	Piercable TPE Capband-8 Yellow in Capmat format	10 mats
1775-3035	Piercable TPE Capband-12 Natural in Capmat format -- please enquire about additional colors	10 mats
	Capclusters	
1775-3001	Piercable TPE Capcluster Natural for capping 96 tubes	10 mats
1775-3004	Piercable TPE Capcluster Blue for capping 96 tubes	10 mats
1775-3005	Piercable TPE Capcluster Green for capping 96 tubes	10 mats
1775-3006	Piercable TPE Capcluster Grey for capping 96 tubes	10 mats
1775-3007	Piercable TPE Capcluster Pink for capping 96 tubes	10 mats
1775-3008	Piercable TPE Capcluster Brown for capping 96 tubes	10 mats
1775-3009	Piercable TPE Capcluster Black for capping 96 tubes	10 mats
1775-3010	Piercable TPE Capcluster Yellow for capping 96 tubes	10 mats
1775-3036	Piercable TPE Capcluster White for capping 96 tubes	10 mats
1775-3038	Piercable TPE Capcluster Purple for capping 96 tubes	10 mats
1775-3039	Piercable TPE Capcluster Light Blue for capping 96 tubes	10 mats
1775-3037	Piercable TPE Capcluster 8 Color capping 96 tubes: Natural, Blue, Green, Grey, Pink, Brown, Black & Yellow	10 mats
1775-3042	Lyo Capclusters	10 mats
1775-3900	Piercable TPE Capcluster Natural for capping 96 Glass tubes	10 mats
	Split Caps	
1775-3027	Split TPE Capband for capping 8 tubes	10 mats
1775-3026	Split TPE Capcluster Natural for capping 96 individual tubes	10 mats
1775-3028	Split TPE Capcluster Blue for capping 96 individual tubes	10 mats
1775-3029	Split TPE Capcluster Green for capping 96 individual tubes	10 mats
1775-3030	Split TPE Capcluster Grey for capping 96 individual tubes	10 mats
1775-3031	Split TPE Capcluster Pink for capping 96 individual tubes	10 mats
1775-3032	Split TPE Capcluster Brown for capping 96 individual tubes	10 mats
1775-3033	Split TPE Capcluster Black for capping 96 individual tubes	10 mats
1775-3034	Split TPE Capcluster Yellow for capping 96 individual tubes	10 mats
	EVA Capping range	Pack Size
1778-2605	EVA Capband break-off for capping 1 tube	1,000 pieces
1772-2722	EVA Capband for capping 8 tubes	1,000 pieces
1773-3000	EVA Capmat for capping 96 tubes	10 mats
1773-3001	EVA Capmat for capping 96 tubes for CoMo Rack-96 only	10 mats
	Accessories	Case Size
1765-5005	Manual Capmatsealer for 1.4ml tubes	1 unit
1765-5006	Manual Capmatsealer for 0.5ml, 0.75ml and 1.4ml tubes	1 unit
1765-5702	Support for 0.5ml tubes	1 unit
1765-5700	Support for 0.75ml tubes	1 unit
1765-5701	Support for 1.1ml tubes	1 unit
1763-5007	Pneumatic Capper	1 unit
1765-4000	Manual Decapper-1	1 unit
1765-4001	Manual Decapper-8	1 unit
1765-5010	Manual Decapper-8 Aluminium	1 unit
1763-5000	Automatic Decapper-96	1 unit