

## NEW! Desktop Single Tube Reader with LCD Display



Camera-based desktop single tube reader features a convenient LCD display area on the front panel that shows the 2D bar code as it reads the tube. An output log records the 2D data for export via a txt file. Requires a Windows XP computer, 2 USB ports and one serial port. Weighs 4 pounds.

- Reads any 2D Data Matrix coded tube in one second
- Small foot print: 6.5" wide x 5" deep x 4.5" tall

1768-5117

Desktop single tube reader



**MICRONIC**

Traceable Sampling Solutions

Please see [www.usascientific.com](http://www.usascientific.com) for more information on Micronic tubes, capping options, scanners, and tube sorters.