


A Perfect Fit: The Micronic TPE Capping Range Provides an Excellent Seal to Greiner Bio-One 96-Well Plates!

If you use Greiner Bio-One plates in your lab, then you should know that the tight, quality seal that is characteristic of the Micronic Capping Range can also protect the samples you have stored in 96-well plates.

- All Micronic TPE and EVA Caps, Capmats, Capclusters, and Capbands fit the Greiner Bio-One 96-well plates (HTS Microplates, Polypropylene plates, Deep Well Storage Plates)
- Improved Micronic TPE Cap material is easier to pierce and better at resealing after piercing
- Micronic TPE Split-Caps are pre-scored and are easily pierced with pipet tips
- TPE Capclusters are available in **11** colors
- Micronic also offers Manual Decappers to assist with cap removal
- All Micronic Caps are standard RNase, DNase, and pyrogen free, and are also available sterile (gamma irradiated)



Micronic Manual Decappers allow you to easily remove a single TPE cap, or a row of 8 TPE caps at a time



Micronic EVA Capmat-96s also provide an excellent seal for Greiner Bio-One 96-well plates



Micronic TPE Capclusters are available in **11** colors for easy sample distinction



The Micronic capping range provides an excellent seal to Micronic racks, as well as to the wide assortment of Greiner Bio-One 96-well

TPE Capping Range Pricing:

Item #	Description	Pack Qty.
Single Caps		
1775-3003	Pierceable TPE Cap for capping 1 tube bulk	960 caps
Capmats		
1775-3000	Pierceable TPE Capmat--96 Natural	10 pcs
1775-3011	Pierceable TPE Capmat--96 Blue	10 pcs
1775-3012	Pierceable TPE Capmat--96 Green	10 pcs
1775-3013	Pierceable TPE Capmat--96 Grey	10 pcs
1775-3014	Pierceable TPE Capmat--96 Pink	10 pcs
1775-3015	Pierceable TPE Capmat--96 Brown	10 pcs
1775-3016	Pierceable TPE Capmat--96 Black	10 pcs
1775-3017	Pierceable TPE Capmat--96 Yellow	10 pcs
Capbands		
1775-3002	Pierceable TPE Capband-8 Natural in Capmat format	10 pcs
1775-3018	Pierceable TPE Capband-8 Blue in Capmat format	10 pcs
1775-3019	Pierceable TPE Capband-8 Green in Capmat format	10 pcs
1775-3020	Pierceable TPE Capband-8 Grey in Capmat format	10 pcs
1775-3021	Pierceable TPE Capband-8 Pink in Capmat format	10 pcs
1775-3022	Pierceable TPE Capband-8 Brown in Capmat format	10 pcs
1775-3023	Pierceable TPE Capband-8 Black in Capmat format	10 pcs
1775-3024	Pierceable TPE Capband-8 Yellow in Capmat format	10 pcs
1775-3035	Pierceable TPE Capband-12 Natural in Capmat format (please inquire about additional Capband-12 colors)	10 pcs
Capclusters		
1775-3001	Pierceable TPE Capcluster-96 Natural	10 pcs
1775-3004	Pierceable TPE Capcluster-96 Blue	10 pcs
1775-3005	Pierceable TPE Capcluster-96 Green	10 pcs
1775-3006	Pierceable TPE Capcluster-96 Grey	10 pcs
1775-3007	Pierceable TPE Capcluster-96 Pink	10 pcs
1775-3008	Pierceable TPE Capcluster-96 Rust	10 pcs
1775-3009	Pierceable TPE Capcluster-96 Black	10 pcs
1775-010	Pierceable TPE Capcluster-96 Yellow	10 pcs
1775-3036	Pierceable TPE Capcluster-96 White	10 pcs
1775-3038	Pierceable TPE Capcluster-96 Purple	10 pcs
1775-3039	Pierceable TPE Capcluster-96 Light Blue	10 pcs
1775-3037	Pierceable TPE Capcluster-96, Case of 8 Colors (Natural, Blue, Green, Grey, Pink, Brown, Black, Yellow)	
1775-3005	Pierceable TPE Capcluster-96, Green (sterile)	24 packs
Split Caps		
1775-3027	Split TPE Capcluster-8	10 pcs
1775-3026	Split TPE Capcluster-96 Natural	10 pcs
1775-3028	Split TPE Capcluster-96 Blue	10 pcs
1775-3029	Split TPE Capcluster-96 Green	10 pcs
1775-3030	Split TPE Capcluster-96 Grey	10 pcs
1775-3031	Split TPE Capcluster-96 Pink	10 pcs
1775-3032	Split TPE Capcluster-96 Brown	10 pcs
1775-3033	Split TPE Capcluster-96 Black	10 pcs
1775-3034	Split TPE Capcluster-96 Yellow	10 pcs
Accessories		
1775-4000	Manual Decapper-1	1 pc
1775-4001	Manual Decapper-8	1 pc
1775-5010	Aluminum Manual Decapper-8	1 unit



Traceable Sampling Solutions

Please see www.usascientific.com for more information on Micronic tubes, capping options, scanners, and tube sorters.