

Sealing Films and Foils

Selection Guide by Attribute

Click Below for Website Link:

	Functional Temperature Range	Product Attribute								Material Construction & Dimensions		
		Sterile	DNase- & RNase- Free	Chemically Resistant	Pierceable	Optically Clear	Transparent	Light Blocking	Breathable	Film / Foil Thickness	Adhesive Thickness	Usable Dimensions
2921-7800 TempPlate® RT Select qPCR Sealing Film	-80°C to +105°C	*	*		*					Polyolefin 50.8µm / 2.0mil	Silicone 50.8µm / 2.0mil	4.574" x 3.083" 116.2mm x 78.3mm
2921-7810 TempPlate® RT Select qPCR Sealing Film, Str	-80°C to +105°C	*	*	*	*					Polyolefin 50.8µm / 2.0mil	Silicone 50.8µm / 2.0mil	4.574" x 3.083" 116.2mm x 78.3mm
2978-2700 TempPlate® RT qPCR Sealing Film, Inset Cut	-40°C to +120°C		*		*					Polyester 50.8µm / 2.0mil	Acrylic 29.2µm / 1.15mil	4.650" x 3.063" 118.1mm x 77.8mm
2978-2100 TempPlate® RT qPCR Sealing Film, Std Cut	-40°C to +120°C		*		*					Polyester 50.8µm / 2.0mil	Acrylic 29.2µm / 1.15mil	4.799" x 3.125" 121.9mm x 79.4mm
2998-1800 TempPlate® EXT Sealing Film Strips	-20°C to +120°C		*				*			Polyester 50.8µm / 2.0mil	Acrylic 41.9µm / 1.65mil	2.885" x 0.396" 73.3mm x 10.1mm
2998-0000 TempPlate® EXT Sealing Film, Inset Cut	-20°C to +120°C		*				*			Polyester 50.8µm / 2.0mil	Acrylic 41.9µm / 1.65mil	4.65" x 3.083" 118.1mm x 78.3mm
2998-4500 TempPlate® EXT Sealing Film, Std Cut	-20°C to +120°C		*				*			Polyester 50.8µm / 2.0mil	Acrylic 41.9µm / 1.65mil	4.80" x 3.208" 121.9mm x 81.5mm
2928-0100 TempPlate® XP Sealing Film, Inset Cut	-40°C to +120°C		*				*			Polyester 50.8µm / 2.0mil	Acrylic 41.9µm / 1.65mil	4.65" x 3.083" 118.1mm x 78.3mm
2972-2100 TempPlate® XP Sealing Film, Standard Cut	-40°C to +120°C		*				*			Polyester 50.8µm / 2.0mil	Acrylic 41.9µm / 1.65mil	4.924" x 3.125" 125.1mm x 79.4mm
2921-0000 TempPlate® Sealing Film	-40°C to +120°C		*				*			Polypropylene 50.8µm / 2.0mil	Acrylic 24.1µm / 0.95mil	4.845" x 3.125" 123.1mm x 79.4mm
2921-0010 TempPlate® Sealing Film, Sterile	-40°C to +120°C	*	*				*			Polypropylene 50.8µm / 2.0mil	Acrylic 24.1µm / 0.95mil	4.845" x 3.125" 123.1mm x 79.4mm
2923-1800 TempPlate® Sealing Foil Strips	-20°C to +120°C		*		*			*		Aluminum foil 50.8µm / 2.0mil	Acrylic 36.8µm / 1.45mil	3.197" x 0.396" 81.2mm x 10.1mm
2923-0100 TempPlate® Sealing Foil	-80°C to +120°C		*	*	*			*		Aluminum foil 35.6µm / 1.4mil	Acrylic 27.9µm / 1.1mil	4.938" x 3.25" 125.4mm x 82.6mm
2923-0110 TempPlate® Sealing Foil, Sterile	-80°C to +120°C	*	*	*	*			*		Aluminum foil 35.6µm / 1.4mil	Acrylic 27.9µm / 1.1mil	4.938" x 3.25" 125.4mm x 82.6mm
2939-6100 Single Tab Sealing Foil, Inset Cut	-80°C to +120°C		*		*			*		Aluminum foil 38.1µm / 1.5mil	Acrylic 37.6µm / 1.48mil	4.625" x 3.063" 117.5mm x 77.8mm
2938-4100 Single Tab Sealing Foil, Std Cut	-80°C to +120°C		*		*			*		Aluminum foil 38.1µm / 1.5mil	Acrylic 37.6µm / 1.48mil	4.875" x 3.25" 123.8mm x 82.6mm
2998-7100 TempPlate® EXT Sealing Foil, Inset Cut	-80°C to +120°C		*	*	*			*		Aluminum foil 48.3µm / 1.9mil	Acrylic 30.5µm / 1.2mil	4.65" x 3.083" 118.1mm x 78.3mm
2998-0100 TempPlate® EXT Sealing Foil, Std Cut	-80°C to +120°C		*	*	*			*		Aluminum foil 48.3µm / 1.9mil	Acrylic 30.5µm / 1.2mil	4.90" x 3.208" 124.5mm x 81.5mm
2920-0000 Microplate Sealing Film	-40°C to +120°C							*		Polyester 50.8µm / 2.0mil	Acrylic 26.7µm / 1.05mil	4.924" x 3.125" 125.1mm x 79.4mm
2920-0010 Microplate Sealing Film, Sterile	-40°C to +120°C	*						*		Polyester 50.8µm / 2.0mil	Acrylic 26.7µm / 1.05mil	4.924" x 3.125" 125.1mm x 79.4mm
2978-5877 EasyApp Microplate Film Kit	-40°C to +120°C							*		Polyester 38.1µm / 1.5mil	Acrylic 26.7µm / 1.05mil	4.75" x 3.125" 120.7mm x 79.4mm
2977-6772 EasyApp Microplate Film Rolls	-40°C to +120°C							*		Polyester 38.1µm / 1.5mil	Acrylic 26.7µm / 1.05mil	4.75" x 3.125" 120.7mm x 79.4mm
2920-3500 Microplate Film, Roll Format	-40°C to +120°C							*		Polyester 50.8µm / 2.0mil	Acrylic 26.7µm / 1.05mil	4.80" x 3.125" 121.9mm x 79.4mm
2920-1010 Microporous Film, Sterile	-20°C to +80°C	*			*				*	Rayon / Acrylic 139.7µm / 5.5mil		5.00" x 3.25" 127.0mm x 82.6mm
2978-5827 Microporous Film Kit, Sterile	-20°C to +80°C	*			*				*	Rayon / Acrylic 139.7µm / 5.5mil		4.75" x 3.25" 120.7mm x 82.6mm
2977-6282 Microporous Film Rolls, Sterile	-20°C to +80°C	*			*				*	Rayon / Acrylic 139.7µm / 5.5mil		4.75" x 3.25" 120.7mm x 82.6mm