

Standard Well Depth



| | | | | | |
|------------------------------------|---------------------------------|--------------------------------------|------------------------------|--------------------------|-------------------|
| Skirt | Semi | Semi | Semi | None | None |
| Additional Features | Contoured sides, rimmed wells | Raised contoured sides, rimmed wells | Straight sides, rimmed wells | Rimmed wells | Chimney wells |
| Notched Corner | A12 | A12 | A12 | A12 | H12 |
| Working Volume (per well) | 1-240 μ l | 1-240 μ l | 1-240 μ l | 1-240 μ l | 1-240 μ l |
| Maximum Volume | 300 μ l | 300 μ l | 300 μ l | 300 μ l | 300 μ l |
| Color Options | Natural, Assorted, White, Black | Natural, Assorted | Natural | Natural, Assorted, White | Natural, Assorted |
| High Contrast Alphanumerics | Yes | Yes | Yes | Yes | Yes |
| Natural Plate # | 1402-9700 | 1402-9300 | 1402-9200 | 1402-9596 | 1402-9600 |
| Catalog Page # | C5 | C5 | C4 | C7 | C6 |

Low Profile Wells



| | | | | | |
|------------------------------------|-------------------------------|-------------------------------------|------------------------------|----------------|-------------------------|
| Skirt | Semi | Semi | Semi | None | Full |
| Additional Features | Contoured sides, rimmed wells | Raised straight sides, rimmed wells | Straight sides, rimmed wells | Rimmed wells | Rimmed wells |
| Notched Corner | A12 | A1 | H12 | H12 | Natural A12 White H1 |
| Working Volume (per well) | 1-140 μ l | 1-140 μ l | 1-140 μ l | 1-140 μ l | 1-140 μ l |
| Maximum Volume | 200 μ l | 200 μ l | 200 μ l | 200 μ l | 200 μ l |
| Color Options | Natural | Natural, White | Natural, White | Natural, White | Natural, White |
| High Contrast Alphanumerics | Yes | Yes | Yes | Yes | Yes |
| Natural Plate # | 1402-9400 | 1402-9100 | 1402-9900 | 1402-9500 | 1402-9800 |
| Catalog Page # | C5 | C4 | C4 | C7 | C8 |

Are you working with very low volumes? Try our low profile PCR tubes or plates.